
pymepix Documentation

Release 1.1.dev0

CFEL Controlled Molecule Imaging group

May 31, 2023

Contents

| | | |
|-----------|--|-----------|
| 1 | Introduction | 3 |
| 2 | Getting started | 5 |
| 2.1 | Installing | 5 |
| 2.1.1 | Installing from PyPI (platform-independent) | 5 |
| 2.1.2 | Installing from git source directly (platform-independent) | 5 |
| 2.1.3 | Build Documentation | 5 |
| 2.2 | Dependencies | 5 |
| 3 | Connecting and Configuring | 7 |
| 3.1 | Connecting | 7 |
| 3.2 | Configuring | 7 |
| 4 | Acquisition | 9 |
| 4.1 | Polling | 9 |
| 4.2 | Callback | 10 |
| 4.3 | Pipelines | 10 |
| 5 | Data Formats | 13 |
| 5.1 | UDP Packets | 13 |
| 5.2 | Decoded Pixels | 13 |
| 5.3 | Decoded Triggers | 14 |
| 5.4 | Time of Flight/Event | 14 |
| 5.5 | Centroid Data | 14 |
| 6 | Examples | 17 |
| 7 | PymepixAcq - Command line | 21 |
| 8 | Pymepix postprocessing | 23 |
| 9 | Troubleshooting | 25 |
| 10 | Pymepix Developer documentation | 27 |
| 10.1 | API reference | 27 |
| 10.1.1 | General overview | 27 |
| 10.1.1.1 | pymepix | 27 |
| 10.1.1.2 | timepixdevice and timepixdef | 27 |

| | | |
|-----------|------------------------------|-----------|
| 10.1.1.3 | config module | 27 |
| 10.1.1.4 | core module | 28 |
| 10.1.1.5 | processing module | 28 |
| 10.1.1.6 | SPIDR module | 28 |
| 10.1.1.7 | util module | 28 |
| 10.1.1.8 | Class overview | 28 |
| 10.1.2 | pymepix | 28 |
| 10.1.2.1 | pymepix package | 28 |
| 11 | References | 67 |
| 12 | Pymepix Documentation | 69 |
| | Bibliography | 71 |
| | Python Module Index | 73 |
| | Index | 75 |

Pymepix documentation

CHAPTER 1

Introduction

Pymepix is intended to bridge the gap between Timepix3 and Python. The goal of the library is to allow a user without deep technical knowledge of Timepix3 to establish a connection, start acquisition, and retrieve and plot pixel and timing information in as few lines of code as possible; at the same time it provides all details to unleash the full power of Timepix3-SPIDR hardware. This is achieved by classes that act as a black-box, handling all of the low level TCP communication and decoding of the UDP data-stream, presenting them in a pythonic fashion. More advanced and lower-level control of SPIDR and Timepix3 is still available from these black-box classes or can be established directly by the user. For easy installation, it only depends on the standard python library, numpy and scikit-learn.

CHAPTER 2

Getting started

2.1 Installing

2.1.1 Installing from PyPI (platform-independent)

Execute `pip install pymepix`. This should install pymepix including all dependencies.

2.1.2 Installing from git source directly (platform-independent)

You can clone pymepix from our main git repository:

```
git clone https://github.com/CFEL-CMI/pymepix.git
```

Navigate into the pymepix library (`cd pymepix`) and run `pip install .`

2.1.3 Build Documentation

To build the documentation for pymepix locally perform the following commands. The first line is only required if there are changes in the package structure or new classes or packages have been added. To only build the existing documentation only the second line must be executed.

```
1 sphinx-apidoc -o ./doc/source/ ./pymepix  
2 python setup.py build_sphinx
```

Adapt `pymepix/config/default.yaml` according to your setup.

2.2 Dependencies

The majority of pymepix only depends on numpy. To use centroiding, the scikit-learn package is required

- *numpy*
- *scikit-learn*: Centroiding and data reduction (Using DBSCAN algorithm for clustering)
- *scipy*: Calculation of the centroids properties from the identified clusters
- *pyzmq*: Inter process communication in the processing pipeline
- *h5py*: Saving processed data as hdf5 files
- *tqdm*: Display a progressbar for post processing
- *pyyaml*: Konfiguration of camera (ip, port, ...)
- *pyserial* (optional): Only used for inclusion of USBTrainID at FLASH and XFEL

CHAPTER 3

Connecting and Configuring

3.1 Connecting

For the camera to work you will have to set up the IP address on your machine, that the camera then communicates with. For Timepix3 with 10 Gb/s that is 192.168.100.1. Look up the official documentation for your camera to find out more.

Before using Pymepix, make sure your camera works properly with the SoPhy software.

The IP address of your TPX camera is the one seen on the OLED screen. Connecting to SPIDR can be done with:

```
>>> timepix = Pymepix(('192.168.100.10', 50000))
```

The number of devices can be found using:

```
>>> len(timepix)  
1
```

Meaning we have one device. To access this device directly, use:

```
tpx0 = timepix[0]
```

And to check the device name:

```
>>> tpx0.deviceName  
W0026_K08
```

3.2 Configuring

To set the biasVoltage to 50 Volts in spidr you can do:

```
>>> timepix.biasVoltage = 50
```

Setting the we can manage its settings directly. To easily setup the device we can use a SoPhy config file (.spx):

```
tpx0.loadConfig('myFile.spx')
```

This sets up all the DAC setting and pixel configurations. Individual parameters can also be set for example. To set the fine threshold to 100 mV do:

```
>>> tpx0.Vthreshold_fine = 100
```

pixel threshold configurations can be set by passing a 256x256 numpy array:

```
import numpy as np  
tpx0.pixelThreshold[...] = 0
```

The same for pixel masks, to set a checkboard mask do:

```
tpx0.pixelMask[::-2] = 1
```

These need to be uploaded to timepix before they take effect:

```
>>> tpx0.uploadPixels()
```

The full list of parameters that can be set can be found in `timepixdevice()`.

CHAPTER 4

Acquisition

Acquisition can be started and stopped by:

```
1 import time
2 from pymepix import Pymepix
3
4 #Connect
5 timepix = Pymepix('192.168.1.10', 50000)
6
7 #Start acquisition
8 timepix.start()
9
10 #Wait
11 time.sleep(1.0)
12
13 #Stop acquisition
14 timepix.stop()
```

Pymepix provides data as a tuple given by (MessageType,data). These are explained in [Data Formats](#). Retrieving the data can be done in two ways: Polling or Callback

4.1 Polling

Polling is where pymepix will place anything retrieved from Timepix into a ring polling buffer. This is the default mode but to reenable it you can use:

```
>>> timepix.enablePolling(maxlen=1000)
```

where *maxlen* describes the maximum number of elements in the buffer before older values are overwritten.

The user can retrieve this data by using:

```
>>> timepix.poll()
(MessageType.RawData, (array[98732405897234589802345, dtype=uint8], 12348798))
```

If there is nothing in the polling buffer then a `PollBufferEmpty` exception is raised. The poll buffer is limited in size but can be extended by doing:

```
>>> timepix.pollBufferLength = 5000
```

This will clear all objects using the polling buffer.

4.2 Callback

The callback method allows the user to deal with the data immediately when it is received. Setting this will clear the polling buffer of any contents.

To set a callback, first you need a function, for example:

```
def my_callback(data_type, data):
    print('My callback is running!!!!')
```

The format of the function must accept two parameters, `MessageType` and an extra data parameter. These are explained in [Data Formats](#). Now to make pymepix use it simply do:

```
>>> timepix.dataCallback = my_callback
```

Now when acquisition is started:

```
>>> timepix.start()
```

The output seen is:

```
... code-block:: sh
```

```
My callback is running!!!! My callback is running!!!! My callback is running!!!! My callback is running!!!! My callback is running!!!!
```

4.3 Pipelines

Pymepix uses pipelines objects in order to process data. Each pipeline is set for each timepix device so each timepix can have a different data pipeline. You can configure them to postprocess or output data in certain ways. For example the `PixelPipeline` object will read from a UDP packet stream and decode the stream into *pixel x*, *pixel y*, *time of arrival* and *time over threshold* arrays. All data is propagated forward through the pipeline so both UDP packets and decoded pixels are output.

To use the (default) `PixelPipeline` pipeline on the first connected timepix device you can do:

```
from pymepix.processing import PixelPipeline, CentroidPipeline

timepix[0].setupAcquisition(PixelPipeline)
```

If you need centroid you instead can do:

```
>>> timepix[0].setupAcquisition(CentroidPipeline)
```

Configuring the pipelines can be done using the `acquisition` property for the timepix device, for example to enable TOFs you can do:

```
>>> timepix[0].acquisition.enableEvents = True
```

A list of pipelines and setting can be found in `acquisition()`

CHAPTER 5

Data Formats

Contains a list of possible data formats output during acquisition. Each entry of the data section represents another element in the tuple. Example shows how to read the data through polling

5.1 UDP Packets

Data Type: MessageType.RawData

Data:

array(uint64) list of UDP packets
uint64 global timer from Timepix at time packets were received

Example:

```
1 data_type, data = timepix.poll()  
2 if data_type is MessageType.RawData:  
3     packets, longtime = data
```

5.2 Decoded Pixels

Data Type: MessageType.PixelData

Data:

array(uint64) pixel x position
array(uint64) pixel y position
array(float) global time of arrival in seconds
array(uint64) time over threshold in nanoseconds

Example:

```
1 data_type, data = timepix.poll()
2 if data_type is MessageType.PixelData:
3     x, y, toa, tot = data
```

5.3 Decoded Triggers

Data Type: MessageType.TriggerData

Data:

- array(uint64)** trigger number
- array(float)** global trigger time in seconds

Example:

```
1 data_type, data = timepix.poll()
2 if data_type is MessageType.TriggerData:
3     t_num, t_time = data
```

5.4 Time of Flight/Event

Data Type: MessageType.EventData

Data:

- array(uint64)** trigger number
- array(uint64)** pixel x position
- array(uint64)** pixel y position
- array(float)** time of flight relative to its trigger in seconds
- array(uint64)** time over threshold in nanoseconds

Example:

```
1 data_type, data = timepix.poll()
2 if data_type is MessageType.EventData:
3     trigger, x, y, tof, tot = data
```

5.5 Centroid Data

Data Type: MessageType.CentroidData

Data:

- array(uint64)** trigger number
- array(float)** center of mass x position
- array(float)** center of mass y position
- array(float)** minimum cluster time of flight
- array(float)** average cluster time over threshold
- array(uint64)** maximum cluster time over threshold

array(uint64) cluster size

Example:

```
1 data_type, data = timepix.poll()
2 if data_type is MessageType.CentroidData:
3     trigger, x, y, tof, avg_tot, max_tot, size = data
```


CHAPTER 6

Examples

Starting timepix and polling data:

```
import pymepix
from pymepix.processing import MessageType
import numpy as np

#Connect to SPIDR
timepix = pymepix.pymepix_connection.PymepixConnection('192.168.1.10', 50000)

#Set bias voltage
timepix.biasVoltage = 50

#Set pixel masks
timepix[0].pixelThreshold = np.zeros(shape=(256,256), dtype=np.uint8)
timepix[0].pixelMask = np.zeros(shape=(256,256), dtype=np.uint8)
timepix[0].uploadPixels()

#Start acquisition
timepix.start()

while True:
    try:
        #Poll
        data_type,data = timepix.poll()
    except pymepix.PollBufferEmpty:
        #If empty then just loop
        continue

    #Handle Raw
    if data_type is MessageType.RawData:

        print('UDP PACKET')

    packets,longtime = data
```

(continues on next page)

(continued from previous page)

```

print('Packet ',packets)
print('Time', longtime)

#Handle Pixels
elif data_type is MessageType.PixelData:

    print('I GOT PIXELS!!!!')

    x,y,toa,tot = data

    print('x',x)
    print('y', y)
    print('toa', toa)
    print('tot',tot)

#Stop
timepix.stop()

```

Using callbacks to acquire:

```

import pymepix
from pymepix.processing import MessageType
import numpy as np
import time

#Connect to SPIDR
timepix = pymepix.Pymepix(('192.168.1.10',50000))

#Set bias voltage
timepix.biasVoltage = 50

#Set pixel masks
timepix[0].pixelThreshold = np.zeros(shape=(256,256),dtype=np.uint8)
timepix[0].pixelMask = np.zeros(shape=(256,256),dtype=np.uint8)
timepix[0].uploadPixels()

#Define callback
def my_callback(data_type,data):
    print('MY CALLBACK!!!!')
    #Handle Raw
    if data_type is MessageType.RawData:

        print('UDP PACKET')

        packets,longtime = data

        print('Packet ',packets)
        print('Time', longtime)

    #Handle Pixels
    elif data_type is MessageType.PixelData:

        print('I GOT PIXELS!!!!')

        x,y,toa,tot = data

```

(continues on next page)

(continued from previous page)

```
print('x',x)
print('y', y)
print('toa', toa)
print('tot',tot)

#Set callback
timepix.dataCallback = my_callback

#Start acquisition
timepix.start()
#Wait 5 seconds
time.sleep(5.0)
#Stop
timepix.stop()
```


CHAPTER 7

PymepixAcq - Command line

Included with pymepix is a command line code using the pymepix library to acquire from timepix. The command line interface

- “connect”: to connect to a running timepix camera and record data
- “post-process”: to post-process recorded raw data files into easier usable hdf5 files containing raw and centroided data

Doing:

```
pymepix-acq --help
```

Outputs the help:

```
usage: pymepix-acq [-h] {connect,post-process} ...

Timepix acquisition script

positional arguments:
  {connect,post-process}
    connect            Connect to TimePix camera and acquire data.
    post-process       Perform post-processing for an acquired raw data file.

optional arguments:
  -h, --help           show this help message and exit
```

You can access the documentation for both commands by executing “pymepix-acq connect -h” or “pymepix-acq post-process -h” respectively.

CHAPTER 8

Pymepix postprocessing

The raw data acquired from the camera could be processed from command line with the command. The processing can also be triggered from the *PymepixViewer*.

Doing:

```
pymepix post-process -f FILE -o OUTPUT_FILE [-t TIMEWALK_FILE] [-c CENT_TIMEWALK_  
FILE] [-n NUMBER_OF_PROCESSES]
```

The generated output file has HDF data format may contain the following datagroups in its root:

- **centroided**
- **raw**
- **timing/timepix**
- **triggers**

The **centroided** datagroups contains the data after centroiding processing. It consists of several datasets : “trigger nr”, “x”, “y”, “tof”, “tot avg”, “tot max”, “clustersize”. Where “trigger nr” is event number, “x”/”y” - coordinates of centroid, tof is time-of-flight (time-of-arrival corrected to the timewalk effect), “tot avg” average value of tot for all voxels in the cluster, “tot max” - max tot value, “clustersize” - the number of voxels in the detected cluster.

The **raw** datagroups contains event data - voxel data with tof synchronized to first trigger. it consists of following datasets: “trigger nr”, “x”, “y”, “tof”, “tot”.

The **timing/timepix** datagroup has only two datasets: “trigger nr”, “timestamp”. Where “trigger nr” contains triggering event numbers from first trigger, while dataset “time” contains the timestamps for the corresponding trigger event in nanosecond in absolute time from the timer of the camera.

Datagroup **triggers** may contain two subgroups “trigger1” and “trigger2” corresponding to the first and second trigger of the camera. Each subgroup consists of only one dataset “time”. These are firing times of the corresponding trigger starting from acquisition in seconds. In case of first trigger these are the times of rising front of the detected trigger pulse. For the second trigger both rising and falling pulse edges are detected. Negative values correspond to the falling edge.

Here's an example to retrieve the data from the HDF5 file into a Pandas DataFrame:

CHAPTER 9

Troubleshooting

- **Whenever there are problems when working with the camera** First make sure you can ping the Timepix camera to ensure a working connection.

Next try starting the SoPhy software and see if it can communicate properly with the camera. Remember to close SoPhy afterwards, as there can only be one process using the address.

- **Make sure to load the correct config file.** If the parameters are off, you might not be able to see anything.
- **Use a flashlight!** When the camera is properly connected and set up, you may use a flashlight to shine directly into the lens and next to it in quick succession.
- **You can see ToA but no ToF data** Check and maybe reconfigure the trigger.
- **Error: Address is already in use** If you get this error, look for any other process that is running and uses the corresponding IP and port. Also try restarting the camera
 - genindex
 - modindex
 - search

CHAPTER 10

Pymepix Developer documentation

This developer documentation contains the API reference for *pymepix*.

10.1 API reference

10.1.1 General overview

The main Pymepix library is built up in several different submodules which each tackle a different task in working with the Timepix camera. As seen in the API index those are * config * core * processing * SPIDR * util

The top layer Pymepix consists of *pymepix*, *timepixdef* and *timepixdevice*.

10.1.1.1 **pymepix**

pymepix provides the highest level of interaction with the library. A single *Pymepix* object will hold all connected Timepix devices and manage the users' interaction with those.

10.1.1.2 **timepixdevice and timepixdef**

A *timepixdevice* object holds all the communication with a single camera. It basically configures, starts and stops. The *timepixdef* has multiple enums to encode all kinds of parameters for Timepix.

10.1.1.3 **config module**

The *config* module gets the information for config parameters.

timepixconf is the base class for the possible configurations.

defaultconfig holds hardcoded config parameters for the camera to initialize.

sophyconfig imports information from SoPhy (.spx) config files. It reads and transforms that information to be used by Pymepix.

10.1.1.4 core module

The *core* module consists of only the log class. It defines functionality for Pymepix' needs and uses the basic python logging module.

10.1.1.5 processing module

The *processing* module provides the data pipeline to process the incoming camera data. Pymepix can use different acquisition pipelines to process the data. Those are defined in *acquisition* with the base functionality provided by *baseacquisition*. An acquisition pipeline determines which steps work in what order on the incoming data and connects those.

Each pipeline consists of acquisition stages (*baseacquisition*), where one stage holds the information about one logical step in the pipeline. Those tasks are currently *udpsampler* (capturing the packets), *rawtodisk* (saving the raw data), *pipeline_packet_processor* (interpreting the raw packets) and *pipeline_centroid_calculator* (compress data by finding blob centers). Each of these specific pipeline steps overwrites the *BasePipelineObject*, which is in fact a python *multiprocessing.Process*.

The majority of the logic for the *pipeline_packet_processor* and the *pipeline_centroid_calculator* is separated in the classes *centroid_calculator* and *packet_processor*. The *pipeline_* classes only add functionality for the integration of those classes into the multiprocessing pipeline.

Each stage knows the task it has to fulfill and then creates one or multiple processes to work on that task in parallel. *datatype*s provides an enum to classify the data that is passed through the pipeline at each step.

10.1.1.6 SPIDR module

This module communicates with the SPIDR chip of the Timepix. One *spidrcontroller* knows about one or more *spidrdevices*. *spidrcmds* lists the known commands to pass information and instructions to a chip. *spidrdefs* extends those commands by constants that can be passed. *error* contains information on possible errors from SPIDR.

10.1.1.7 util module

storage provides some functionality to save data.

spidrDummyTCP and *-UDP* can be used to simulate a timepix camera. Both are still rudimentary but helpful for debugging.

spidrDummyTCP accepts packets in so the configuration of timepix can be tested.

spidrDummyUDP samples and sends packets from a given file into the void. This can be used to test the pipeline functionality by capturing those packets with Pymepix.

10.1.1.8 Class overview

10.1.2 pymepix

10.1.2.1 pymepix package

Subpackages

pymepix.SPIDR package

Submodules

pymepix.SPIDR.error module

```
exception pymepix.SPIDR.error.PymePixException(error_code)
Bases: Exception

ERR_STR = ['no error', 'ERR_UNKNOWN_CMD', 'ERR_MSG_LENGTH', 'ERR_SEQUENCE', 'ERR_ILLEGAL']
MONITOR_ERR_STR = ['MON_ERR_TEMP_DAQ', 'MON_ERR_POWER_DAQ']
SPIDR_ERR_STR = ['SPIDR_ERR_I2C_INIT', 'SPIDR_ERR_LINK_INIT', 'SPIDR_ERR_MPL_INIT', 'SPIDR_ERR_SC_INIT']
STORE_ERR_STR = ['no error', 'STORE_ERR_TPX', 'STORE_ERR_WRITE', 'STORE_ERR_WRITE_CHECK']
TPX3_ERR_STR = ['no error', 'TPX3_ERR_SC_ILLEGAL', 'TPX3_ERR_SC_STATE', 'TPX3_ERR_SC_ERROR']
errorMessage(code)

class pymepix.SPIDR.error.SPIDRErrorDfs
Bases: object

ERR_ADC_HARDW = 7
ERR_DAC_HARDW = 8
ERR_FLASH_STORAGE = 10
ERR_ILLEGAL_PAR = 4
ERR_MONITOR = 11
ERR_MON_HARDW = 9
ERR_MSG_LENGTH = 2
ERR_NONE = 0
ERR_NOT_IMPLEMENTED = 5
ERR_SEQUENCE = 3
ERR_TPX3_HARDW = 6
ERR_UNKNOWN_CMD = 1
```

pymepix.SPIDR.spidrcmds module

This module contains a list of all (found) commands for the SPIDR board

```
class pymepix.SPIDR.spidrcmds.SpidrCmds
Bases: enum.IntEnum

A class that packages all the commands under a single name

CMD_AUTOTRIG_START = 1090
CMD_AUTOTRIG_STOP = 1091
CMD_BIAS_SUPPLY_ENA = 1375
CMD_BURN_EFUSE = 297
CMD_CLEAR_BUSY = 2313
CMD_CONFIG_CTPR = 288
```

```
CMD_DDRIVEN_READOUT = 1094
CMD_DECODERS_ENA = 1377
CMD_DISPLAY_INFO = 2315
CMD_ERASE_ADDRPORTS = 1652
CMD_ERASE_DAC = 1653
CMD_ERASE_PIXCONF = 1655
CMD_ERASE_REGISTERS = 1654
CMD_GET_ADC = 1352
CMD_GET_AVDD = 1355
CMD_GET_AVDD_NOW = 1357
CMD_GET_BOARDID = 2318
CMD_GET_CHIPBOARDID = 2319
CMD_GET_CTPR = 290
CMD_GET_DAC = 282
CMD_GET_DEVICECOUNT = 2317
CMD_GET_DEVICEID = 272
CMD_GET_DEVICEIDS = 273
CMD_GET_DEVICEPORT = 278
CMD_GET_DVDD = 1356
CMD_GET_DVDD_NOW = 1359
CMD_GET_EFUSES = 296
CMD_GET_EXTSHUTTERCNTR = 1366
CMD_GET_FANSPEED = 1385
CMD_GET_FIRMWAREVERSION = 2306
CMD_GET_FPGATEMP = 1384
CMD_GET_GENCONFIG = 820
CMD_GET_GPIO = 1920
CMD_GET_HEADERFILTER = 2309
CMD_GET_HUMIDITY = 1390
CMD_GET_IPADDR_DEST = 276
CMD_GET_IPADDR_SRC = 274
CMD_GET_LOCALTEMP = 1354
CMD_GET_OUTBLOCKCONFIG = 828
CMD_GET_PIXCONF = 557
CMD_GET_PLLCONFIG = 822
CMD_GET_PRESSURE = 1391
```

```
CMD_GET_PWRPULSECONFIG = 1371
CMD_GET_READOUTSPEED = 1380
CMD_GET_REMOTETEMP = 1353
CMD_GET_SERVERPORT = 279
CMD_GET_SHUTTERCNTR = 1367
CMD_GET_SHUTTEREND = 1365
CMD_GET_SHUTTERSTART = 1364
CMD_GET_SLVSCONFIG = 830
CMD_GET_SOFTWVERSION = 2305
CMD_GET_SPIDRREG = 1923
CMD_GET_SPIDR_ADC = 1358
CMD_GET_STARTOPTS = 1661
CMD_GET_TIMER = 1362
CMD_GET_TPNUMBER = 819
CMD_GET_TPPERIODPHASE = 816
CMD_GET_TRIGCONFIG = 1088
CMD_GET_VDD = 1388
CMD_GET_VDD_NOW = 1389
CMD_MASK = 65535
CMD_NOP = 0
CMD_NOREPLY = 524288
CMD_PAUSE_READOUT = 1095
CMD_PWRPULSE_ENA = 1373
CMD_READ_FLASH = 1662
CMD_REINIT_DEVICE = 294
CMD_REINIT_DEVICES = 295
CMD_REPLY = 65536
CMD_RESET_COUNTERS = 1368
CMD_RESET_DEVICE = 292
CMD_RESET_DEVICES = 293
CMD_RESET_MODULE = 2311
CMD_RESET_PIXELS = 558
CMD_RESET_TIMER = 1361
CMD_RESTART_TIMERS = 1360
CMD_SELECT_CHIPBOARD = 1387
CMD_SEQ_READOUT = 1093
```

```
CMD_SET_BIAS_ADJUST = 1376
CMD_SET_BOARDID = 1926
CMD_SET_BUSY = 2312
CMD_SET_CHIPBOARDID = 1925
CMD_SET_CTPR = 289
CMD_SET_CTPR_LEON = 291
CMD_SET_DAC = 283
CMD_SET_DACS_DFLT = 287
CMD_SET_EXTDAC = 826
CMD_SET_FANSPEED = 1386
CMD_SET_GENCONFIG = 821
CMD_SET_GPIO = 1921
CMD_SET_GPIO_PIN = 1922
CMD_SET_HEADERFILTER = 2310
CMD_SET_IPADDR_DEST = 277
CMD_SET_IPADDR_SRC = 275
CMD_SET_LOGLEVEL = 2314
CMD_SET_OUTBLOCKCONFIG = 829
CMD_SET_OUTPUTMASK = 1378
CMD_SET_PIXCONF = 554
CMD_SET_PLLCONFIG = 823
CMD_SET_PWRPULSECONFIG = 1372
CMD_SET_READOUTSPEED = 1379
CMD_SET_SENSEDAC = 824
CMD_SET_SERVERPORT = 281
CMD_SET_SLVSCONFIG = 831
CMD_SET_SPIDRREG = 1924
CMD_SET_TIMEOFDAY = 2316
CMD_SET_TIMER = 1363
CMD_SET_TPNUMBER = 818
CMD_SET_TPPERIODPHASE = 817
CMD_SET_TRIGCONFIG = 1089
CMD_STORE_ADDRPORTS = 1648
CMD_STORE_DACS = 1649
CMD_STORE_PIXCONF = 1651
CMD_STORE_REGISTERS = 1650
```

```
CMD_STORE_STARTOPTS = 1660
CMD_T0_SYNC = 1381
CMD_TPX_POWER_ENA = 1374
CMD_UPLOAD_PACKET = 827
CMD_VALID_ADDRPORTS = 1656
CMD_VALID_DACS = 1657
CMD_VALID_PIXCONF = 1659
CMD_VALID_REGISTERS = 1658
CMD_WRITE_FLASH = 1663
```

pymepix.SPIDR.spidrcontroller module

SPIDR related classes

class `pymepix.SPIDR.spidrcontroller.SPIDRController(dst_ip_port, src_ip_port)`
Bases: `pymepix.core.log.Logger`

Object that interfaces over ethernet with the SPIDR board

This object interfaces with the spidr board through TCP and is used to send commands and receive data. It can be treated as a list of SpidrDevice objects to talk to a specific device

Parameters

- `dst_ip_port` (tuple of str and int) – socket style tuple of SPIDR ip address and port
- `src_ip_port` (tuple of str and int, optional) – socket style tuple of the IP address and port of the interface that is connecting to SPIDR

Examples

The class can be used to talk to SPIDR

```
>>> spidr = SPIDRController(('192.168.1.10', 50000))
>>> spidr.fpgaTemperature
39.5
```

Or access a specific SpidrDevice (e.g. Timepix/Medipix)

```
>>> spidr[0].deviceId
7272
>>> spidr[1].deviceId
2147483648
```

Warning: This object assumes SPIDR is working as intended however since this is still in development there are a few functions that do not behave as they should, this will be documented in their relevant areas.

CpuToTpx

Cpu2Tpx register access

Parameters `value` (`int`) – Value to write to the register
Returns Current value of the register
Return type int
Raises PymePixException – Communication error

Notes

Register controls clock setup

DeviceAndPorts

ShutterTriggerCount

Number of times the shutter is triggered in auto trigger mode

Parameters `value` (`int`) – Trigger count to set for auto trigger mode (Set to 0 for infinite triggers)
Returns Current value of the trigger count read from SPIDR
Return type int
Raises PymePixException – Communication error

ShutterTriggerCtrl

Shutter Trigger Control register access

Parameters `value` (`int`) – Value to write to the register
Returns Current value of the register
Return type int
Raises PymePixException – Communication error

ShutterTriggerDelay

Delay time before shutter can be triggered again in auto trigger mode

Parameters `value` (`int`) – Time in ns
Returns value – Current time in ns read from SPIDR
Return type int
Raises PymePixException – Communication error

ShutterTriggerFreq

Triggering frequency for the auto trigger

Parameters `value` (`float`) – Frequency in mHz
Returns Frequency value in mHz read from SPIDR
Return type float
Raises PymePixException – Communication error

ShutterTriggerLength

Length of time shutter remains open at each trigger

Parameters `value` (`int`) – Length in ns
Returns value – Current length in ns read from SPIDR
Return type int

Raises PymePixException – Communication error

ShutterTriggerMode

Controls how the shutter is triggered

Parameters **value** (SpidrShutterMode) – Shutter trigger mode to set

Returns Current shutter operation mode read from SPIDR

Return type SpidrShutterMode

Raises PymePixException – Communication error

Notes

AutoTrigger is the only functioning trigger mode that SPIDR can operate in

TdcTriggerCounter

Trigger packets sent by SPIDR since last counter reset

UdpMonPacketCounter

UdpPacketCounter

UDP packets sent by SPIDR since last counter reset

UdpPausePacketCounter

UDP packets collected during readout pause since last counter reset

avdd

avddNow

biasVoltage

Bias voltage

Parameters **volt** (*int*) – Bias voltage to supply in volts Minimum is 12V and Maximum is 104V

Returns Current bias supply in volts

Return type int

Raises PymePixException – Communication error

chipboardFanSpeed

chipboardId

clearBusy ()

closeShutter ()

Immediately closes the shutter

Raises PymePixException – Communication error

convertHtonl (*x*)

convertNtohl (*x*)

dataDrivenReadout ()

Set SPIDR into data driven readout mode

Data driven mode refers to the pixels packets sent as they are hit rather than camera style frames

Raises PymePixException – Communication error

Warning: This is the only tested mode for pymepix. It is recommended that this is enabled

deviceCount

Count of devices connected to SPIDR

Returns Number of devices connected to SPIDR

Return type int

Raises PymePixException – Communication error

Warning: SPIDR always returns 4 since it currently can't determine if the devices are actually valid or not

deviceIds

The ids of all devices connected to the SPIDR board

Returns A list all connected device ids

Return type list of int

Raises PymePixException – Communication error

Notes

Index of devices are the same as the those in the SPIDRController list

```
>>> spidr[1].deviceId == spidr.deviceIds[1]
True
```

disableExternalRefClock()

SPIDR receives its reference clock internally

This should be set in single SPIDR mode. When combining other SPIDR board, the master will set this to disabled

Raises PymePixException – Communication error

disablePeriphClk80Mhz()

dvdd

dvddNow

enableDecoders(*enable*)

Determines whether the internal FPGA decodes ToA values

Time of Arrival from UDP packets are gray encoded if this is enabled then SPIDR will decode them for you, otherwise you have to do this yourself after extracting them

Parameters **enable** (bool) – True - enable FPGA decoding False - disable FPGA decoding

Raises PymePixException – Communication error

Tip: Enable this

enableExternalRefClock()

SPIDR receives its reference clock externally

This is often used when combining multiple Timepixs together so they can synchronize their clocks. The SPIDR board essentially acts as a slave to other SPIDRs

Raises PymePixException – Communication error

enablePeriphClk80Mhz()**externalShutterCounter****firmwareVersion**

Firmware version

Returns Version number of firmware within the FPGA

Return type int

Raises PymePixException – Communication error

fpgaTemperature

Temperature of FPGA board read from sensor

Returns Temperature in Celsius

Return type float

Raises PymePixException – Communication error

getAdc(channel, nr_of_samples)**getSpidrReg(addr)****humidity**

Humidity read from sensor

Returns Humidity as percentage

Return type int

Raises PymePixException – Communication error

linkCounts**localTemperature**

Local ?????!?!? Temperature read from sensor

Returns Temperature in Celsius

Return type float

Raises PymePixException – Communication error

openShutter()

Immediately opens the shutter indefinitely

Raises PymePixException – Communication error

Notes

This overwrites shutter configurations with one that forces an open shutter

pauseReadout()**pressure**

Pressure read from sensor

Returns Pressure in bar

Return type int

Raises PymePixException – Communication error

reinitDevices()

Resets and initializes all devices

Raises PymePixException – Communication error

remoteTemperature

Remote ????!!? Temperature read from sensor

Returns Temperature in Celsius

Return type float

Raises PymePixException – Communication error

request(cmd, dev_nr, message_length, expected_bytes=0)

Sends a command and (may) receive a reply

Parameters

- **cmd** (SpidrCmds) – Command to send
- **dev_nr** (int) – Device to send the request to. 0 is SPIDR and device number n is n+1
- **message_length** (int) – Length of the message in bytes
- **expected_bytes** (int) – Length of expected reply from request (if any) (Default: 0)

Returns Returns a numpy array of ints if reply expected, otherwise None

Return type numpy.array of int or None

Raises PymePixException – Communication error

requestGetBytes(cmd, dev_nr, expected_bytes, args=0)

requestGetInt(cmd, dev_nr, arg=0)

requestGetIntBytes(cmd, dev_nr, expected_bytes, args=0)

requestGetInts(cmd, dev_nr, num_ints, args=0)

requestSetInt(cmd, dev_nr, value)

requestSetIntBytes(cmd, dev_nr, value_int, value_bytes)

requestSetInts(cmd, dev_nr, value)

resetCounters()

resetDevices()

Resets all devices

resetModule(readout_speed)

Resets the SPIDR board and sets a new readout speed

Parameters **readout_speed** (SpidrReadoutSpeed) – Read-out speed the device will operate at

Notes

Its not clear if this does anything as its not usually used

resetPacketCounters ()

resetTimers ()

Resets all timers to zero

Sets the internal 48-bit timers for all Timepix/Medipix devices to zero

Raises PymePixException – Communication error

restartTimers ()

Restarts SPIDR and Device timers

Synchronizes both the SPIDR clock and Timepix/Medipix clocks so both trigger and ToA timestamps match

Important: This must be done if event selection is required (e.g. time of flight) otherwise the timestamps will be offset

Raises PymePixException – Communication error

sequentialReadout (tokens, now)

setBiasSupplyEnable (enable)

Enables/Disables bias supply voltage

Parameters **enable** (*bool*) – True - enables bias supply voltage False - disables bias supply voltage

Raises PymePixException – Communication error

setBusy ()

setPowerPulseEnable (enable)

setShutterTriggerConfig (mode, length_us, freq_hz, count, delay_ns=0)

Set the shutter configuration in one go

Parameters

- **mode** (*int*) – Shutter trigger mode
- **length_us** (*int*) – Shutter open time in microseconds
- **freq_hz** (*int*) – Auto trigger frequency in Hertz
- **count** (*int*) – Number of triggers
- **delay_ns** (*int, optional*) – Delay between each trigger (Default: 0)

Raises PymePixException – Communication error

setSpidrReg (addr, value)

setTpxPowerPulseEnable (enable)

shutterCounter

shutterTriggerConfig

```
softwareVersion
    Software version

        Returns Version number of software in the SPIDR board

        Return type int

        Raises PymePixException – Communication error

spidrFanSpeed

startAutoTrigger()
    Starts the auto trigger

        Raises PymePixException – Communication error

stopAutoTrigger()
    Stops the auto trigger

        Raises PymePixException – Communication error

vdd

vddNow

pymepix.SPIDR.spidrcontroller.main()
```

pymepix.SPIDR.spidrdefs module

Module that contains constants that can be passed into spidr

```
class pymepix.SPIDR.spidrdefs.SpidrReadoutSpeed
    Bases: enum.Enum

    An enumeration.

    Default = 0
    HighSpeed = 2309737967
    LowSpeed = 305419896

class pymepix.SPIDR.spidrdefs.SpidrRegs
    Bases: enum.IntEnum

    An enumeration.

    SPIDR_CPU2TPX_WR_I = 456
    SPIDR_DEVICES_AND_PORTS_I = 704
    SPIDR_FE_GTX_CTRL_STAT_I = 768
    SPIDR_IPMUX_CONFIG_I = 896
    SPIDR_PIXEL_FILTER_I = 916
    SPIDR_PIXEL_PKTOUNTER_I = 832
    SPIDR_PIXEL_PKTOUNTER_OLD_I = 912
    SPIDR_SHUTTERTRIG_CNT_I = 660
    SPIDR_SHUTTERTRIG_CTRL_I = 656
    SPIDR_SHUTTERTRIG_DELAY_I = 684
```

```

SPIDR_SHUTTERTRIG_FREQ_I = 664
SPIDR_SHUTTERTRIG_LENGTH_I = 668
SPIDR_TDC_TRIGGERCOUNTER_I = 760
SPIDR_UDPMON_PKTNUMBER_I = 904
SPIDR_UDPPAUSE_PKTNUMBER_I = 908
SPIDR_UDP_PKTNUMBER_I = 900

class pymepix.SPIDR.spidrdefs.SpidrShutterMode
Bases: enum.Enum

An enumeration.

Auto = 4
ExternalFallingRising = 1
ExternalFallingTimer = 3
ExternalRisingFalling = 0
ExternalRisingTimer = 2
Open = 6
PulseCounter = 5

```

pymepix.SPIDR.spidrdevice module

```

class pymepix.SPIDR.spidrdevice.SpidrDevice(spidr_ctrl, device_num)
Bases: pymepix.core.log.Logger

Object that interfaces with a specific device (Timepix/Medipix) connect to SPIDR

This object handles communication and management of a specific device. There is no need to create this object directly as SpidrController automatically creates it for you and is accessed by its [] getter methods

```

Parameters

- **spidr_ctrl** (SpidrController) – SPIDR controller object the device belongs to
- **device_num** – Device index from SPIDR (Starts from 1)

TpPeriodPhase

clearPixelConfig()

columnTestPulseRegister

deviceId

Returns unique device Id

Parameters

- **spidr_ctrl** (SpidrController) – SPIDR controller object the device belongs to
- **device_num** – Device index from SPIDR (Starts from 1)

devicePort

genConfig

getDac(dac_code)

```
getDacOut (nr_samples)
getPixelConfig ()
headerFilter
ipAddrDest
ipAddrSrc
linkStatus
outBlockConfig
pixelPacketCounter
pllConfig
powerPulseConfig
readoutSpeed
reinitDevice ()
reset ()
resetPixelConfig (index=-1, all_pixels=False)
resetPixels ()
serverPort
setDac (dac_code, dac_val)
setDacDefault ()
setExternalDac (dac_code, dac_val)
setHeaderFilter (eth_mask, cpu_mask)
setOutputMask (value)
setPixelMask (mask)
setPixelTestBit (test)
setPixelThreshold (threshold)
setSenseDac (dac_code)
setSinglePixelMask (x, y, mask)
setSinglePixelTestBit (x, y, val)
setSinglePixelThreshold (x, y, threshold)
setTpPeriodPhase (period, phase)
shutterEnd
shutterStart
slaveConfig
t0Sync ()
timer
tpNumber
uploadPacket (packet)
```

```
pymepix.uploadPixelConfig(formatted=True, columns_per_packet=1)
```

Module contents

pymepix.clustering package

Submodules

pymepix.clustering.cluster_stream module

```
class pymepix.clustering.cluster_stream.ClusterStream(dim=256, max_dist_tof=1e-08, min_cluster_size=3, tot_offset=0.5, *args, **kwargs)
```

Bases: object

```
perform(data)
```

Module contents

pymepix.config package

Submodules

pymepix.config.defaultconfig module

```
class pymepix.config.defaultconfig.DefaultConfig
Bases: pymepix.config.timepixconfig.TimepixConfig
```

Provides default values for DAC parameters

```
biasVoltage()
```

Returns bias Voltage

```
dacCodes()
```

Accessor for the dac parameters

Returns The value for every DAC parameter

Return type list of tuples (<dac code>, <value>)

```
maskPixels
```

Returns mask pixels

```
testPixels
```

Returns test pixels

```
thresholdPixels
```

Returns threshold pixels

pymepix.config.load_config module

```
pymepix.config.load_config.load_config(config_name='default.yaml')
```

pymepix.config.sophyconfig module

```
class pymepix.config.sophyconfig.SophyConfig(filename)
    Bases: pymepix.config.timepixconfig.TimepixConfig, pymepix.core.log.Logger

This class provides functionality for interpreting a .spx config file from SoPhy.

biasVoltage()
    Returns bias Voltage

dacCodes()
    Accessor for the dac parameters

        Returns The value for every DAC parameter

        Return type list of tuples (<dac code>, <value>)

filename

loadFile(filename)

maskPixels
    Accessor for the mask pixels [0, 1]

        Returns The information which pixels are to be masked

        Return type numpy.ndarray (256, 256)

parseDAC(xmlstring)
    Reads and formats DAC parameters

parsePixelConfig(zip_file, file_names)
    Reads and formats the pixel data from config file.
```

Notes

The spx config file saves the pixel information row by row while the timepix camera expects the information column wise.

```
saveMask()

testPixels
    Accessor for the test pixels

        Returns

        Return type numpy.ndarray (256, 256)

thresholdPixels
    Accessor for the pixel thresholds [0, 15]

        Returns The threshold information for each pixel

        Return type numpy.ndarray (256, 256)
```

pymepix.config.sophyconfig.main()

pymepix.config.timepixconfig module

```
class pymepix.config.timepixconfig.TimepixConfig
    Bases: abc.ABC
```

```
biasVoltage()
    Returns bias Voltage

dacCodes()
    Returns an iterator with format daccode,value

maskPixels
    Returns mask pixels

testPixels
    Returns test pixels

thresholdPixels
    Returns threshold pixels
```

Module contents

pymepix.core package

Submodules

pymepix.core.log module

```
class pymepix.core.log.Logger(name)
    Bases: pymepix.core.log.PymepixLogger
    Standard logging using logger library

        Parameters name (str) – Name used for logging

    getLogger(name)

class pymepix.core.log.ProcessLogger(name)
    Bases: pymepix.core.log.PymepixLogger
    Sends logs to queue to be processed by logging thread

        Parameters name (str) – Name used for logging

    getLogger(name)
```

Module contents

pymepix.processing package

Subpackages

pymepix.processing.logic package

Submodules

pymepix.processing.logic.centroid_calculator module

```
class pymepix.processing.logic.centroid_calculator.CentroidCalculator(cent_timewalk_lut=None,
                                                                     num-
                                                                     ber_of_processes=4,
                                                                     clus-
                                                                     ter-
                                                                     ing_args={},
                                                                     db-
                                                                     scan_clustering=True,
                                                                     *args,
                                                                     **kwargs)
```

Bases: *pymepix.processing.logic.processing_step.ProcessingStep*

Class responsible for calculating centroids in timepix data. This includes the calculation of the clusters first and the centroids. The data processed is not the direct raw data but the data that has been processed by the PacketProcessor before (x, y, tof, tot).

process(data) :

Process data and return the result. To use this class only this method should be used! Use the other methods only for testing or if you are sure about what you are doing

calculate_centroids_cluster_stream(chunk)

calculate_centroids_dbscan(chunk)

calculate_centroids_properties(shot, x, y, tof, tot, labels)

Calculates the properties of the centroids from labeled data points.

ATTENTION! The order of the points can have an impact on the result due to errors in the floating point arithmetics.

Very simple example: arr = np.random.random(100) arr.sum() - np.sort(arr).sum() This example shows that there is a very small difference between the two sums. The inaccuracy of floating point arithmetics can depend on the order of the values. Strongly simplified (3.2 + 3.4) + 2.7 and 3.2 + (3.4 + 2.7) can be unequal for floating point numbers.

Therefore there is no guarantee for strictly equal results. Even after sorting. The error we observed can be about 10^-22 nano seconds.

Currently this issue exists only for the TOF-column as the other columns are integer-based values.

centroid_chunks_to_centroids(chunks)

centroids = [[] for i in range(7)] for chunk in list(chunks):

if chunk != None:

for index, coordinate in enumerate(chunk): centroids[index].append(coordinate)

cluster_stream_preprocess(shot, x, y, tof, tot)

cs_max_dist_tof

Setting the maximal ToF distance between the voxels belonging to the cluster in Cluster Streaming algorithm

cs_min_cluster_size

Setting the minimal cluster size in Cluster Streaming algorithm

cs_sensor_size

Setting for the number of packets skipped during processing. Every packet_skip packet is processed. This means for a value of 1 every packet is processed. For 2 only every 2nd packet is processed.

cs_tot_offset

Setting the ToT ratio factor of the voxel to the ToT of previous voxel in Cluster Streaming algorithm. Zero factor means ToT of prev. voxel should be larger. 0.5 factor means ToT of prev voxel could be high than the half of the considered voxel

dbscan_clustering**epsilon****min_samples****perform_centroiding_cluster_stream**(*chunks*)**perform_centroiding_dbSCAN**(*chunks*)**perform_clustering_dbSCAN**(*shot, x, y, tof*)

The clustering with DBSCAN, which is performed in this function is dependent on the order of the data in rare cases. Therefore, reordering in any means can lead to slightly changed results, which should not be an issue.

Martin Ester, Hans-Peter Kriegel, Jiří Sander, Xiaowei Xu: A Density Based Algorithm for Discovering Clusters [p. 229-230] (<https://www.aaai.org/Papers/KDD/1996/KDD96-037.pdf>) A more specific explanation can be found here: <https://stats.stackexchange.com/questions/306829/why-is-dbscan-deterministic>

process(*data*)**tot_threshold**

Determines which time over threshold values to filter before centroiding

This is useful in reducing the computational time in centroiding and can filter out noise.

triggers_processed

Setting for the number of packets skipped during processing. Every packet_skip packet is processed. This means for a value of 1 every packet is processed. For 2 only every 2nd packet is processed.

```
class pymepix.processing.logic.centroid_calculator.CentroidCalculatorPooled(number_of_processes,  
                           *args,  
                           **kwargs)
```

Bases: *pymepix.processing.logic.centroid_calculator.CentroidCalculator*

Parallelized implementation of CentroidCalculator using mp.Pool for parallelization.

perform_centroiding(*chunks*)**post_process**()**pre_process**()

```
pymepix.processing.logic.centroid_calculator.calculate_centroids_dbSCAN(chunk,  
                           tot_threshold,  
                           _tof_scale,  
                           ep-  
                           silon,  
                           min_samples,  
                           _cent_timewalk_lut)
```

```
pymepix.processing.logic.centroid_calculator.calculate_centroids_properties(shot,  
                                x,  
                                y,  
                                tof,  
                                tot,  
                                la-  
                                bels,  
                                _cent_timewalk_lut)
```

Calculates the properties of the centroids from labeled data points.

ATTENTION! The order of the points can have an impact on the result due to errors in the floating point arithmetics.

Very simple example: arr = np.random.random(100) arr.sum() - np.sort(arr).sum() This example shows that there is a very small difference between the two sums. The inaccuracy of floating point arithmetics can depend on the order of the values. Strongly simplified $(3.2 + 3.4) + 2.7$ and $3.2 + (3.4 + 2.7)$ can be unequal for floating point numbers.

Therefore there is no guarantee for strictly equal results. Even after sorting. The error we observed can be about 10^{-22} nano seconds.

Currently this issue exists only for the TOF-column as the other columns are integer-based values.

```
pymepix.processing.logic.centroid_calculator.perform_clustering_dbSCAN(shot,  
                                x, y,  
                                tof,  
                                _tof_scale,  
                                ep-  
                                silon,  
                                min_samples)
```

The clustering with DBSCAN, which is performed in this function is dependent on the order of the data in rare cases. Therefore, reordering in any means can lead to slightly changed results, which should not be an issue.

Martin Ester, Hans-Peter Kriegel, Jiří Sander, Xiaowei Xu: A Density Based Algorithm for Discovering Clusters [p. 229-230] (<https://www.aaai.org/Papers/KDD/1996/KDD96-037.pdf>) A more specific explanation can be found here: <https://stats.stackexchange.com/questions/306829/why-is-dbscan-deterministic>

pymepix.processing.logic.packet_processor module

```
class pymepix.processing.logic.packet_processor.PacketProcessor(handle_events=True,  
                           event_window=(0.0,  
                                         10000.0), posi-  
                                         tion_offset=(0,  
                                                       0), orienta-  
                                         tion=<PixelOrientation.Up:  
                                         0>,  
                                         start_time=0,  
                                         time-  
                                         walk_lut=None,  
                                         *args,  
                                         **kwargs)
```

Bases: *pymepix.processing.logic.processing_step.ProcessingStep*

Class responsible to transform the raw data coming from the timepix directly into an easier processible data format. Takes into account the pixel- and trigger data to calculate toa and tof dimensions.

process(data) :

Process data and return the result. To use this class only this method should be used! Use the other methods

only for testing or if you are sure about what you are doing

```
clearBuffers()
correct_global_time(arr, ltime)
event_window
find_events_fast()
find_events_fast_post()
    Call this function at the very end of to also have the last two trigger events processed
getBuffers(val_filter=None)
handle_events
    Type noindex
orientPixels(col, row)
    Orient the pixels based on Timepix orientation
post_process()
pre_process()
process(data)
process_pixels(pixdata, longtime)
process_trigger1(pixdata, longtime)
process_trigger2(tidtrigdata, longtime)
updateBuffers(val_filter)
```

class pymepix.processing.logic.packet_processor.**PixelOrientation**
Bases: enum.IntEnum

Defines how row and col are interpreted in the output

```
Down = 2
    x=-column, y = -row
Left = 1
    x=row, y=-column
Right = 3
    x=-row, y=column
Up = 0
    Up is the default, x=column,y=row
```

pymepix.processing.logic.processing_step module

```
class pymepix.processing.logic.processing_step.ProcessingStep(name, parameter_wrapper_class=<class 'pymepix.processing.logic.shared_processor'
```

Bases: *pymepix.core.log.Logger*, abc.ABC

Representation of one processing step in the pipeline for processing timepix raw data. Implementations are provided by PacketProcessor and CentroidCalculator. To combine those (and possibly other) classes into a pipeline they have to implement this interface. Also provides pre- and post-process implementations which are required for integration in the online processing pipeline (see PipelineCentroidCalculator and PipelinePacketProcessor).

Currently the picture is the following:

- For post processing the CentroidCalculator and the PacketProcessor are used directly
- PipelineCentroidCalculator and PipelinePacketProcessor build on top of CentroidCalculator and PacketProcessor to provide an integration in the existing online processing pipeline for online analysis.

```
post_process()  
pre_process()  
process(data)
```

pymepix.processing.logic.shared_processing_parameter module

```
class pymepix.processing.logic.shared_processing_parameter.SharedProcessingParameter(value)  
Bases: object
```

Variang of the ProcessingParameter used for sharing among multiple processes. This class has to be used if running with the multiprocessing pipeline to ensure all instances of the processing classes are updated when parameters are changed.

value

```
exception pymepix.processing.logic.shared_processing_parameter.UnknownParameterTypeException  
Bases: Exception
```

Module contents

Submodules

pymepix.processing.acquisition module

Module that contains predefined acquisition pipelines for the user to use

```
class pymepix.processing.acquisition.CentroidPipeline(data_queue, address, long-  
time)  
Bases: pymepix.processing.acquisition.PixelPipeline
```

A Pixel pipeline that includes centroiding

Same as the pixel pipeline but also includes centroid processing, note that this can be extremely slow when dealing with a huge number of objects

numBlobProcesses

Number of python processes to spawn for centroiding

Setting this will spawn the appropriate number of processes to perform centroiding. Changes take effect on next acquisition.

```
class pymepix.processing.acquisition.PixelPipeline(data_queue, address, longtime,  
use_event=False, name='Pixel',  
event_window=(0, 0.001))  
Bases: pymepix.processing.baseacquisition.AcquisitionPipeline
```

An acquisition pipeline that includes the udpsampler and pixel processor

A pipeline that will read from a UDP address and decode the pixels a useable form. This class can be used as a base for all acquisition pipelines.

pymepix.processing.baseacquisition module

Module deals with managing processing objects to form a data pipeline

class `pymepix.processing.baseacquisition.AcquisitionPipeline(name, data_queue)`

Bases: `pymepix.core.log.Logger`

Class that manages various stages

addStage (`stage_number, pipeline_klass, *args, num_processes=1, **kwargs`)

Adds a stage to the pipeline

getStage (`stage_number`)

isRunning

stages

start()

Starts all stages

stop()

Stops all stages

class `pymepix.processing.baseacquisition.AcquisitionStage(stage,`

`num_processes=1)`

Bases: `pymepix.core.log.Logger`

Defines a single acquisition stage

Usually not created directly. Instead created by `AcquisitionPipeline` Represent a single pipeline stage and handles management of queues and message passing as well as creation and destruction of processing objects.

Processes are not created until build() is called and do not run until start() is called

Parameters `stage (int)` – Initial position in the pipeline, lower stages are executed first

build (`input_queue=None, output_queue=None, file_writer=None`)

configureStage (`pipeline_klass, *args, **kwargs`)

Configures the stage with a particular processing class

Parameters

- `pipeline_klass` (`BasePipeline`) – A pipeline class object
- `*args` – positional arguments to pass into the class init
- `**kwargs` – keyword arguments to pass into the class init

numProcess

Number of processes to spawn when built

Parameters `value (int)` – Number of processes to spawn when acquisition starts

Returns Number of processes

Return type int

outputQueue

processes

setArgs (`*args, **kwargs`)

```
stage
    Current position in the pipeline

start()
startTrainID()
stop(force=False)
stopTrainID()

pymepix.processing.baseacquisition.main()
```

pymepix.processing.basepipeline module

Base implementation of objects relating to the processing pipeline

```
class pymepix.processing.basepipeline.BasePipelineObject(name,           in-
                                                               put_queue=None,
                                                               create_output=True,
                                                               num_outputs=1,
                                                               shared_output=None,
                                                               propogate_input=True)
Bases: multiprocessing.context.Process, pymepix.core.log.ProcessLogger
```

Base class for integration in a processing pipeline

Parameters

- **name** (*str*) – Name used for logging
- **input_queue** (*multiprocessing.Queue*, optional) – Data queue to perform work on (usually) from previous step in processing pipeline
- **create_output** (*bool*, *optional*) – Whether this creates its own output queue to pass data, ignored if (Default: True)
- **num_outputs** (*int*, *optional*) – Used with create_output, number of output queues to create (Default: 1)
- **shared_output** (*multiprocessing.Queue*, optional) – Data queue to pass results into, useful when multiple processes can put data into the same queue (such as results from centroiding). Ignored if create_output is True (Default: None)
- **propogate_input** (*bool*) – Whether the input data should be propgated further down the chain

enable

Enables processing

Determines whether the class will perform processing, this has the result of signalling the process to terminate. If there are objects ahead of it then they will stop receiving data if an input queue is required then it will get from the queue before checking processing This is done to prevent the queue from growing when a process behind it is still working

Parameters **value** (*bool*) – Enable value

Returns Whether the process is enabled or not

Return type *bool*

classmethod hasOutput()

Defines whether this class can output results or not, e.g. Centroiding can output results but file writing classes do not

Returns Whether results are generated

Return type bool

outputQueues

Exposes the outputs so they may be connected to the next step

Returns All of the outputs

Return type list of multiprocessing.Queue

post_run()

Function called after main processing loop, override to

pre_run()

Function called before main processing loop, override to

process (data_type=None, data=None)

Main processing function, override this do perform work

To perform work within the pipeline, a class must override this function. General guidelines include, check for correct data type, and must return None for both if no output is given.

pushOutput (data_type, data)

Pushes results to output queue (if available)

Parameters

- **data_type** (int) – Identifier for data type (see MessageType for types)
- **data** (any) – Results from processing (must be pickleable)

run()

Method to be run in sub-process; can be overridden in sub-class

pymepix.processing.basepipeline.main()

pymepix.processing.datatypes module

Defines data that is passed between processing objects

class pymepix.processing.datatypes.MessageType

Bases: enum.IntEnum

Defines the type of message that is being passed into a multiprocessing queue

CentroidData = 3

Centroided Data

CloseFileCommand = 5

Close File Message

EventData = 2

Event Data

OpenFileCommand = 4

Open File message

PixelData = 1

Decoded Pixel/Trigger Data

```
RawData = 0
    Raw UDP packets

TriggerData = 8
    Decoded Triggers
```

pymepix.processing.pipeline_centroid_calculator module

Processors relating to centroiding

```
class pymepix.processing.pipeline_centroid_calculator.PipelineCentroidCalculator(centroid_calculator):
    pymepix.processing.pipeline_basepipeline.BasePipelineObject
    = <pymepix.processing.pipeline_basepipeline_basepipeline>, in-
    put_queue=None, create_output=True, num_outputs=1, shared_output=None
```

Bases: *pymepix.processing.basepipeline.BasePipelineObject*

Performs centroiding on EventData received from Packet processor

```
process(data_type=None, data=None)
    Main processing function, override this to perform work
```

To perform work within the pipeline, a class must override this function. General guidelines include, check for correct data type, and must return None for both if no output is given.

pymepix.processing.pipeline_packet_processor module

```
class pymepix.processing.pipeline_packet_processor.PipelinePacketProcessor(packet_processor):
    pymepix.processing.log
    = <pymepix.processing.log>, in-
    put_queue=None, create_output=True, num_outputs=1, shared_output=None
```

Bases: *pymepix.processing.basepipeline.BasePipelineObject*

Processes Pixel packets for ToA, ToT, triggers and events

This class, creates a UDP socket connection to SPIDR and receives the UDP packets from Timepix. It then pre-processes them and sends them off for more processing

```
init_new_process()
```

```
post_run()
```

Function called after main processing loop, override to

pre_run()
Function called before main processing loop, override to

process (data_type=None, data=None)
Main processing function, override this do perform work

To perform work within the pipeline, a class must override this function. General guidelines include, check for correct data type, and must return None for both if no output is given.

pymepix.processing.rawfilesampler module

```
class pymepix.processing.rawfilesampler.RawFileSampler(file_name, output_file, number_of_processes=None, timewalk_file=None, cent_timewalk_file=None, progress_callback=None, clustering_args={}, db_scan_clustering=True, **kwargs)
```

Bases: object

bytes_from_file (chunksize=8192)

handle_lsb_time (pixdata)

handle_msb_time (pixdata)

handle_other (pixdata)
trash data which arrives before 1st timestamp data (heartbeat)

init_new_process (file)
create connections and initialize variables in new process

post_run ()

pre_run ()
init stuff which should only be available in new process

push_data (post=False)

run ()
method which is executed in new process via multiprocessing.Process.start

saveToHDF5 (output_file, raw, clusters, timeStamps, trigger1, trigger2)

pymepix.processing.rawtodisk module

```
class pymepix.processing.rawtodisk.Raw2Disk(context=None)
```

Bases: *pymepix.core.log.ProcessLogger*

Class for asynchronously writing raw files Intended to allow writing of raw data while minimizing impact on UDP reception reliability.

close (socket)
Close the file currently in progress. call in main below

open_file (socket, filename)
Creates a file with a given filename and path.

this doesn't work anylonger using 2 sockets for the communication functionality needs to be put outside where you have access to the socket

write (*data*)

Writes data to the file. Parameter is buffer type (e.g. bytearray or memoryview)

Not sure how useful this function actually is... It completes the interface for this class but from a performance point of view it doesn't improve things. How could this be benchmarked?

`pymepix.processing.rawtodisk.main()`

`pymepix.processing.rawtodisk.main_process()`

seperate process not strictly necessary, just to double check if this also works with multiprocessing doesn't work for debugging

pymepix.processing.udpsampler module

class `pymepix.processing.udpsampler.UdpSampler` (*address*, *longtime*, *chunk_size*=10000, *flush_timeout*=0.3, *input_queue*=None, *create_output*=True, *num_outputs*=1, *shared_output*=None)

Bases: `multiprocessing.context.Process`, `pymepix.core.log.ProcessLogger`

Recieves udp packets from SPIDR

This class, creates a UDP socket connection to SPIDR and recivies the UDP packets from Timepix It them pre-processes them and sends them off for more processing

close_file

create_socket_connection (*address*)

Establishes a UDP connection to spidr

enable

Enables processing

Determines whether the class will perform processing, this has the result of signalling the process to terminate. If there are objects ahead of it then they will stop receiving data if an input queue is required then it will get from the queue before checking processing This is done to prevent the queue from growing when a process behind it is still working

Parameters **value** (*bool*) – Enable value

Returns Whether the process is enabled or not

Return type *bool*

get_useful_packets (*packet*)

init_new_process ()

create connections and initialize variables in new process

outfile_name

post_run ()

method get's called either at the very end of the process live or if there's a socket timeout and raw2disk file should be closed

pre_run ()

init stuff which should only be available in new process

record

Enables saving data to disk

Determines whether the class will perform processing, this has the result of signalling the process to terminate. If there are objects ahead of it then they will stop receiving data if an input queue is required then it will get from the queue before checking processing. This is done to prevent the queue from growing when a process behind it is still working

Parameters `value (bool)` – Enable value

Returns Whether the process should record and write to disk or not

Return type bool

run()

method which is executed in new process via multiprocessing.Process.start

pymepix.processing.udpsampler.**main()**

pymepix.processing.usbtrainid module

class pymepix.processing.usbtrainid.**USBTrainID** (`name='USBTrainId'`)

Bases: multiprocessing.context.Process, *pymepix.core.log.ProcessLogger*

Class for asynchronously writing raw files. Intended to allow writing of raw data while minimizing impact on UDP reception reliability

connect_device(device)

Establish connection to USB device

run()

Method to be run in sub-process; can be overridden in sub-class

pymepix.processing.usbtrainid.**main()**

Module contents

pymepix.util package

Submodules

pymepix.util.spidrDummyTCP module

class pymepix.util.spidrDummyTCP.**TPX3Handler** (`request, client_address, server`)

Bases: socketserver.BaseRequestHandler, *pymepix.core.log.Logger*

handle()

pymepix.util.spidrDummyTCP.**main()**

pymepix.util.spidrDummyUDP module

pymepix.util.spidrDummyUDP.**main()**

pymepix.util.storage module

Useful functions to store data

```
pymepix.util.storage.open_output_file(filename, ext, index=0)
pymepix.util.storage.store_centroid(f, data)
pymepix.util.storage.store_raw(f, data)
pymepix.util.storage.store_toa(f, data)
pymepix.util.storage.store_tof(f, data)
```

pymepix.util.tcpsampler module

```
class pymepix.util.tcpsampler.TcpSampler(address,      longtime,      chunk_size=10000,
                                            flush_timeout=0.3,      input_queue=None,
                                            create_output=True,      num_outputs=1,
                                            shared_output=None)
Bases: multiprocessing.context.Process, pymepix.core.log.ProcessLogger
```

Recieves tcp packets

The same as UdpSampler just with TCP

close_file

createConnection(address)

Establishes a TCP connection

enable

Enables processing

Determines whether the class will perform processing, this has the result of signalling the process to terminate. If there are objects ahead of it then they will stop receiving data if an input queue is required then it will get from the queue before checking processing This is done to prevent the queue from growing when a process behind it is still working

Parameters value (bool) – Enable value

Returns Whether the process is enabled or not

Return type bool

get_useful_packets(packet)

init_new_process()

outfile_name

post_run()

pre_run()

record

Enables saving data to disk

Determines whether the class will perform processing, this has the result of signalling the process to terminate. If there are objects ahead of it then they will stop receiving data if an input queue is required then it will get from the queue before checking processing This is done to prevent the queue from growing when a process behind it is still working

Parameters value (bool) – Enable value

Returns Whether the process should record and write to disk or not

Return type bool

```
run()
    Method to be run in sub-process; can be overridden in sub-class

stopRaw2Disk()
    self.debug('Stopping Raw2Disk') self.write2disk.close() self.write2disk.my_sock.send_string('SHUTDOWN')
    # print(write2disk.my_sock.recv()) self.write2disk.write_thr.join() self.debug('Raw2Disk stopped')

pymepix.util.tcpsampler.main()
```

pymepix.util.timewalk module

```
pymepix.util.timewalk.compute_timewalk(tof, tot, region)
pymepix.util.timewalk.compute_timewalk_lookup(tof, tot, region)
```

Module contents

Submodules

pymepix.main module

Main module for pymepix

```
pymepix.main.connect_timepix(args)
pymepix.main.main()
pymepix.main.post_process(args)
```

pymepix.post_processing module

```
class pymepix.post_processing.ProgressBar(iterator=None, desc=None, total=None,
                                         leave=True, file=None, ncols=None, mininterval=0.1, maxinterval=10.0, miniters=None,
                                         ascii=None, disable=False, unit='it',
                                         unit_scale=False, dynamic_ncols=False,
                                         smoothing=0.3, bar_format=None, initial=0, position=None, postfix=None,
                                         unit_divisor=1000, write_bytes=False,
                                         lock_args=None, nrows=None, colour=None,
                                         delay=0, gui=False, **kwargs)
```

Bases: tqdm.std.tqdm

```
gui_bar_fun = None
update_to(progress)
```

```
pymepix.post_processing.run_post_processing(input_file_name, output_file,
                                              number_processes, time-
                                              walk_file, cent_timewalk_file,
                                              progress_callback=<function updatePro-
                                              gressBar>, clustering_args={}, db-
                                              scan_clustering=True, **kwargs)
pymepix.post_processing.updateProgressBar(progress)
```

pymepix.pymepix_connection module

exception pymepix.pymepix_connection.PollBufferEmpty

Bases: Exception

```
class pymepix.pymepix_connection.PymepixConnection(spidr_address=('192.168.1.10',
                                                               50000),
                                                       src_ip_port=('192.168.1.1',
                                                               8192), pipeline_class=<class
'pymepix.processing.acquisition.PixelPipeline'>)
```

Bases: *pymepix.core.log.Logger*

High level class to work with timepix and perform acquisition

This class performs connection to SPIDR, initialization of timepix and handling of acquisition. Each individual timepix device can be accessed using the square bracket operator.

Parameters

- **spidr_address** (tuple of str and int) – socket style tuple of SPIDR ip address and port
- **src_ip_port** (tuple of str and int, optional) – socket style tuple of the IP address and port of the interface that is connecting to SPIDR

Examples

Startup device

```
>>> timepix = Pymepix(('192.168.1.10', 50000))
```

Find how many Timepix are connected

```
>>> len(timepix)
1
```

Set the Bias voltage >>> timepix.biasVoltage = 50

Access a specific Timepix device:

```
>>> timepix[0].deviceName
W0026_K06
```

Load a config file into timepix

```
>>> timepix[0].loadSophyConfig('W0026_K06_50V.spx')
```

biasVoltage

Bias voltage in volts

dataCallback

Function to call when data is received from a timepix device

This has the effect of disabling polling.

data_thread()**enablePolling (maxlen=100)**

Enables polling mode

This clears any user defined callbacks and the polling buffer

getDevice (num) → pymepix.timepixdevice.TimepixDevice**isAcquiring****numDevices****poll (block=False)**

If polling is used, returns data stored in data buffer.

the buffer is in the form of a ring and will overwrite older values if it becomes full

Returns

Return type MessageType , data

pollBufferLength

Get/Set polling buffer length

Clears buffer on set

start()

Starts acquisition

stop()

Stops acquisition

pymepix.timepixdef module

```
class pymepix.timepixdef.DacRegisterCodes
Bases: enum.IntEnum
```

An enumeration.

Ibias_CP_PLL = 17

Ibias_DiscS1_OFF = 9

Ibias_DiscS1_ON = 8

Ibias_DiscS2_OFF = 11

Ibias_DiscS2_ON = 10

Ibias_Ikrum = 4

Ibias_PixelDAC = 12

Ibias_Preamp_OFF = 2

Ibias_Preamp_ON = 1

Ibias_TPbufferIn = 13

Ibias_TPbufferOut = 14

```
PLL_Vcntrl = 18
VPreamp_NCAS = 3
VTP_coarse = 15
VTP_fine = 16
Vfbk = 5
Vthreshold_coarse = 7
Vthreshold_fine = 6

class pymepix.timepixdef.GrayCounter
Bases: enum.IntEnum
An enumeration.

Disable = 0
Enable = 8
Mask = 8

class pymepix.timepixdef.OperationMode
Bases: enum.IntEnum
An enumeration.

EventiTot = 4
Mask = 6
ToA = 2
ToAandToT = 0

class pymepix.timepixdef.PacketType
Bases: enum.Enum
An enumeration.

Pixel = 1
Trigger = 0

class pymepix.timepixdef.Polarity
Bases: enum.IntEnum
An enumeration.

Negative = 1
Positive = 0

class pymepix.timepixdef.SuperPixel
Bases: enum.IntEnum
An enumeration.

Disable = 0
Enable = 64
Mask = 64
```

```
class pymepix.timepixdef.TestPulse
Bases: enum.IntEnum

An enumeration.

Disable = 0
Enable = 32
Mask = 32

class pymepix.timepixdef.TestPulseDigAnalog
Bases: enum.IntEnum

An enumeration.

DiscriminatorDigital = 512
FrontEndAnalog = 0
Mask = 512

class pymepix.timepixdef.TestPulseGenerator
Bases: enum.IntEnum

An enumeration.

External = 1024
Internal = 0
Mask = 1024

class pymepix.timepixdef.TimeofArrivalClock
Bases: enum.IntEnum

An enumeration.

Mask = 2048
PhaseShiftedGray = 0
SystemClock = 2048

class pymepix.timepixdef.TimerOverflow
Bases: enum.IntEnum

An enumeration.

CycleOverflow = 0
Mask = 128
StopOverflow = 128
```

pymepix.timepixdevice module

```
exception pymepix.timepixdevice.BadPixelFormat
Bases: Exception

exception pymepix.timepixdevice.ConfigClassException
Bases: Exception
```

```
class pymepix.timepixdevice.TimepixDevice(spidr_device, data_queue,
                                             pipeline_class=<class
                                             'pymepix.processing.acquisition.PixelPipeline'>)
```

Bases: *pymepix.core.log.Logger*

Provides high level control of a timepix/medipix object

Ibias_DiscS1_OFF

[0, 15]

Ibias_DiscS1_ON

[0, 255]

Ibias_DiscS2_OFF

[0, 15]

Ibias_DiscS2_ON

[0, 255]

Ibias_Ikrum

[0, 255]

Ibias_PixelDAC

[0, 255]

Ibias_Preamp_OFF

[0, 15]

Ibias_Preamp_ON

[0, 255]

Ibias_TPbufferIn

[0, 255]

Ibias_TPbufferOut

[0, 255]

VPreamp_NCAS

[0, 255]

VTP_coarse

[0, 255]

VTP_fine

[0, 511]

Vfbk

[0, 255]

Vthreshold_coarse

[0, 15]

Vthreshold_fine

[0, 511]

acquisition

Returns the acquisition object

Can be used to set parameters in the acquisition directly for example, to setup TOF calculation when using a PixelPipeline

```
>>> tpx.acquisition.enableEvents  
False  
>>> tpx.acquisition.enableEvents = True
```

config

devIdToString()
Converts device ID into readable string

Returns Device string identifier

Return type str

deviceName

grayCounter

loadConfig(*args, **kwargs)
Loads dac settings from the Config class

operationMode

pauseHeartbeat()

pixelMask
Pixel mask set for timepix device

Parameters **value** (numpy.array of int) – 256x256 uint8 threshold mask to set locally

Returns Locally stored pixel mask matrix

Return type numpy.array of int or None

pixelTest
Pixel test set for timepix device

Parameters **value** (numpy.array of int) – 256x256 uint8 pixel test to set locally

Returns Locally stored pixel test matrix

Return type numpy.array of int or None

pixelThreshold
Threshold set for timepix device

Parameters **value** (numpy.array of int) – 256x256 uint8 threshold to set locally

Returns Locally stored threshold matrix

Return type numpy.array of int or None

polarity

refreshPixels()
Loads timepix pixel configuration to local array

resetPixels()
Clears pixel configuration

resumeHeartbeat()

setConfigClass(klass: pymepix.config.timepixconfig.TimepixConfig)

setDac(code, value)
Sets the DAC parameter using codes

Parameters

- **code** (int) – DAC code to set
- **value** (int) – value to set

```
setEthernetFilter(eth_filter)
    Sets the packet filter, usually set to 0xFFFF to all all packets

setupAcquisition(acquisition_klass, *args, **kwargs)

setupDevice()
    Sets up valid paramters for acquisition

    This will be manual when other acquistion parameters are working

start()
start_recording(path)
stop()
stop_recording()

superPixel

testPulse

testPulseDigitalAnalog
testPulseGeneratorSource
timeOfArrivalClock
timerOverflowControl

update_timer()
    Heartbeat thread

uploadPixels()
    Uploads local pixel configuration to timepix

pymepix.timepixdevice.main()
```

Module contents

CHAPTER 11

References

CHAPTER 12

Pymepix Documentation

[Pymepix](#) is a python library for interfacing, controlling and acquiring from SPIDR-Timepix detectors.

See also the accompanying [Pymepix-viewer](#) for an example user tool.

See [\[AlRefaie2019\]](#) for a description of version 1.0 of both tools and as a formal, citeable reference and the [online manual](#) for further information on later versions.

Bibliography

[AlRefaie2019] A. Al-Refaie, M. Johny, J. Correa, D. Pennicard, P. Svihra, A. Nomerotski, S. Trippel, J. Küpper: PymePix: a python library for SPIDR readout of Timepix3, *Journal of Instrumentation* **14**, P10003–P10003 (2019), DOI: [10.1088/1748-0221/14/10/p10003](https://doi.org/10.1088/1748-0221/14/10/p10003); arXiv:1905.07999 [physics]

Python Module Index

p

pymepix, 66
pymepix.clustering, 43
pymepix.clustering.cluster_stream, 43
pymepix.config, 45
pymepix.config.defaultconfig, 43
pymepix.config.load_config, 43
pymepix.config.sophyconfig, 44
pymepix.config.timepixconfig, 44
pymepix.core, 45
pymepix.core.log, 45
pymepix.main, 59
pymepix.post_processing, 59
pymepix.processing, 57
pymepix.processing.acquisition, 50
pymepix.processing.baseacquisition, 51
pymepix.processing.basepipeline, 52
pymepix.processing.datatypes, 53
pymepix.processing.logic, 50
pymepix.processing.logic.centroid_calculator,
 46
pymepix.processing.logic.packet_processor,
 48
pymepix.processing.logic.processing_step,
 49
pymepix.processing.logic.shared_processing_parameter,
 50
pymepix.processing.pipeline_centroid_calculator,
 54
pymepix.processing.pipeline_packet_processor,
 54
pymepix.processing.rawfilesampler, 55
pymepix.processing.rawtodisk, 55
pymepix.processing.udpsampler, 56
pymepix.processing.usbtrainid, 57
pymepix.pymepix_connection, 60
pymepix.SPIDR, 43
pymepix.SPIDR.error, 29
pymepix.SPIDR.spidrcmds, 29
pymepix.SPIDR.spidrcontroller, 33
pymepix.SPIDR.spidrdefs, 40
pymepix.SPIDR.spidrdevice, 41
pymepix.timepixdef, 61
pymepix.timepixdevice, 63
pymepix.util, 59
pymepix.util.spidrDummyTCP, 57
pymepix.util.spidrDummyUDP, 57
pymepix.util.storage, 58
pymepix.util.tcpsampler, 58
pymepix.util.timewalk, 59

Index

A

acquisition (pymepix.timepixdevice.TimepixDevice attribute), 64

AcquisitionPipeline (class in pymepix.processing.baseacquisition), 51

AcquisitionStage (class in pymepix.processing.baseacquisition), 51

addStage () (pymepix.processing.baseacquisition.AcquisitionPipeline method), 51

Auto (pymepix.SPIDR.spidrdefs.SpidrShutterMode attribute), 41

avdd (pymepix.SPIDR.spidrcontroller.SPIDRController attribute), 35

avddNow (pymepix.SPIDR.spidrcontroller.SPIDRController attribute), 35

B

BadPixelFormat, 63

BasePipelineObject (class in pymepix.processing.basepipeline), 52

biasVoltage (pymepix.pymepix_connection.PymepixConnection attribute), 60

biasVoltage (pymepix.SPIDR.spidrcontroller.SPIDRController attribute), 35

biasVoltage () (pymepix.config.defaultconfig.DefaultConfig method), 43

biasVoltage () (pymepix.config.sophyconfig.SophyConfig method), 44

biasVoltage () (pymepix.config.timepixconfig.TimepixConfig method), 44

build () (pymepix.processing.baseacquisition.AcquisitionStage method), 51

bytes_from_file () (pymepix.processing.rawfilesampler.RawFileSampler method), 55

calculate_centroids_dbSCAN () (in module pymepix.processing.logic.centroid_calculator), 47

calculate_centroids_dbSCAN () (pymepix.processing.logic.centroid_calculator.CentroidCalculator method), 46

calculate_centroids_properties () (in module pymepix.processing.logic.centroid_calculator), 47

calculate_centroids_properties () (pymepix.processing.logic.centroid_calculator.CentroidCalculator method), 46

Centroid_chunks_to_centroids () (pymepix.processing.logic.centroid_calculator.CentroidCalculator method), 46

CentroidCalculator (class in pymepix.processing.logic.centroid_calculator), 46

CentroidCalculatorPooled (class in pymepix.processing.logic.centroid_calculator), 47

CentroidData (pymepix.processing.datatypes.MessageType attribute), 53

CentroidPipeline (class in pymepix.processing.acquisition), 50

chipboardFanSpeed (pymepix.SPIDR.spidrcontroller.SPIDRController attribute), 35

chipboardId (pymepix.SPIDR.spidrcontroller.SPIDRController attribute), 35

clearBuffers () (pymepix.processing.logic.packet_processor.PacketProcessor method), 49

clearBusy () (pymepix.SPIDR.spidrcontroller.SPIDRController method), 35

clearPixelConfig () (pymepix.SPIDR.spidrdevice.SpidrDevice

C

C

```

close_file (pymepix.processing.udpsampler.UdpSampler@MD_GET_CHIPBOARDID
           attribute), 56
close_file (pymepix.util.tcpsampler.TcpSampler at-
           tribute), 58
CloseFileCommand (pymepix.processing.datatypes.MessageType attribute), 30
attribute), 53
closeShutter () (pymepix.SPIDR.spidrcontroller.SPIDRController@attribute), 30
method), 35
cluster_stream_preprocess () (pymepix.processing.logic.centroid_calculator.CentroidCalculator@attribute), 30
method), 46
ClusterStream (class in
pymepix.clustering.cluster_stream), 43
CMD_AUTOTRIG_START (pymepix.SPIDR.spidrcmds.SpidrCmds attribute), 29
CMD_AUTOTRIG_STOP (pymepix.SPIDR.spidrcmds.SpidrCmds attribute), 29
CMD_BIAS_SUPPLY_ENA (pymepix.SPIDR.spidrcmds.SpidrCmds attribute), 29
CMD_BURN_EFUSE (pymepix.SPIDR.spidrcmds.SpidrCmds@CMD_GET_EFUSES (pymepix.SPIDR.spidrcmds.SpidrCmds
attribute), 29
CMD_CLEAR_BUSY (pymepix.SPIDR.spidrcmds.SpidrCmds@CMD_GET_EXTSHUTTERCNTR
attribute), 29
CMD_CONFIG_CTPR (pymepix.SPIDR.spidrcmds.SpidrCmds attribute), 29
CMD_DDRIVEN_READOUT (pymepix.SPIDR.spidrcmds.SpidrCmds attribute), 29
CMD_DECODERS_ENA (pymepix.SPIDR.spidrcmds.SpidrCmds attribute), 30
CMD_DISPLAY_INFO (pymepix.SPIDR.spidrcmds.SpidrCmds attribute), 30
CMD_ERASE_ADDRPORTS (pymepix.SPIDR.spidrcmds.SpidrCmds attribute), 30
CMD_ERASE_DACS (pymepix.SPIDR.spidrcmds.SpidrCmds attribute), 30
CMD_ERASE_PIXCONF (pymepix.SPIDR.spidrcmds.SpidrCmds attribute), 30
CMD_ERASE_REGISTERS (pymepix.SPIDR.spidrcmds.SpidrCmds attribute), 30
CMD_GET_ADC (pymepix.SPIDR.spidrcmds.SpidrCmds attribute), 30
CMD_GET_AVDD (pymepix.SPIDR.spidrcmds.SpidrCmds attribute), 30
CMD_GET_AVDD_NOW (pymepix.SPIDR.spidrcmds.SpidrCmds attribute), 30
CMD_GET_BOARDID (pymepix.SPIDR.spidrcmds.SpidrCmds attribute), 30
CMD_GET_CHIPBOARDID (pymepix.SPIDR.spidrcmds.SpidrCmds attribute), 30
CMD_GET_CTPR (pymepix.SPIDR.spidrcmds.SpidrCmds attribute), 30
CMD_GET_DAC (pymepix.SPIDR.spidrcmds.SpidrCmds attribute), 30
CMD_GET_DEVICECOUNT (pymepix.SPIDR.spidrcmds.SpidrCmds attribute), 30
CMD_GET_DEVICEID (pymepix.SPIDR.spidrcmds.SpidrCmds attribute), 30
CMD_GET_DEVICEIDS (pymepix.SPIDR.spidrcmds.SpidrCmds attribute), 30
CMD_GET_DEVICEPORT (pymepix.SPIDR.spidrcmds.SpidrCmds attribute), 30
CMD_GET_DVDD (pymepix.SPIDR.spidrcmds.SpidrCmds attribute), 30
CMD_GET_DVDD_NOW (pymepix.SPIDR.spidrcmds.SpidrCmds attribute), 30
CMD_GET_EFUSES (pymepix.SPIDR.spidrcmds.SpidrCmds attribute), 30
CMD_GET_EXTSHUTTERCNTR (pymepix.SPIDR.spidrcmds.SpidrCmds attribute), 30
CMD_GET_FANSPED (pymepix.SPIDR.spidrcmds.SpidrCmds attribute), 30
CMD_GET_FIRMWAREVERSION (pymepix.SPIDR.spidrcmds.SpidrCmds attribute), 30
CMD_GET_GENCONFIG (pymepix.SPIDR.spidrcmds.SpidrCmds attribute), 30
CMD_GET_FPGATEMP (pymepix.SPIDR.spidrcmds.SpidrCmds attribute), 30
CMD_GET_GPIO (pymepix.SPIDR.spidrcmds.SpidrCmds attribute), 30
CMD_GET_HEADERFILTER (pymepix.SPIDR.spidrcmds.SpidrCmds attribute), 30
CMD_GET_HUMIDITY (pymepix.SPIDR.spidrcmds.SpidrCmds attribute), 30
CMD_GET_IPADDR_DEST (pymepix.SPIDR.spidrcmds.SpidrCmds attribute), 30
CMD_GET_IPADDR_SRC (pymepix.SPIDR.spidrcmds.SpidrCmds attribute), 30
CMD_GET_LOCALTEMP (pymepix.SPIDR.spidrcmds.SpidrCmds attribute), 30
CMD_GET_OUTBLOCKCONFIG (pymepix.SPIDR.spidrcmds.SpidrCmds attribute), 30

```

| | | |
|--|-----|--|
| CMD_GET_PIXCONF (pymepix.SPIDR.spidrcmds.SpidrCmds attribute), 30 | at- | CMD_GET_VDD (pymepix.SPIDR.spidrcmds.SpidrCmds attribute), 31 |
| CMD_GET_PLLCONFIG (pymepix.SPIDR.spidrcmds.SpidrCmds attribute), 30 | at- | CMD_MASK (pymepix.SPIDR.spidrcmds.SpidrCmds attribute), 31 |
| CMD_GET_PRESSURE (pymepix.SPIDR.spidrcmds.SpidrCmds attribute), 30 | at- | CMD_NOP (pymepix.SPIDR.spidrcmds.SpidrCmds attribute), 31 |
| CMD_GET_PWRPULSECONFIG (pymepix.SPIDR.spidrcmds.SpidrCmds attribute), 30 | at- | CMD_NOREPLY (pymepix.SPIDR.spidrcmds.SpidrCmds attribute), 31 |
| CMD_GET_READOUTSPEED (pymepix.SPIDR.spidrcmds.SpidrCmds attribute), 31 | at- | CMD_PAUSE_READOUT (pymepix.SPIDR.spidrcmds.SpidrCmds attribute), 31 |
| CMD_GET_REMOTETEMP (pymepix.SPIDR.spidrcmds.SpidrCmds attribute), 31 | at- | CMD_PWRPULSE_ENA (pymepix.SPIDR.spidrcmds.SpidrCmds attribute), 31 |
| CMD_GET_SERVERPORT (pymepix.SPIDR.spidrcmds.SpidrCmds attribute), 31 | at- | CMD_READ_FLASH (pymepix.SPIDR.spidrcmds.SpidrCmds attribute), 31 |
| CMD_GET_SHUTTERCNTR (pymepix.SPIDR.spidrcmds.SpidrCmds attribute), 31 | at- | CMD_REINIT_DEVICE (pymepix.SPIDR.spidrcmds.SpidrCmds attribute), 31 |
| CMD_GET_SHUTTEREND (pymepix.SPIDR.spidrcmds.SpidrCmds attribute), 31 | at- | CMD_REINIT_DEVICES (pymepix.SPIDR.spidrcmds.SpidrCmds attribute), 31 |
| CMD_GET_SHUTTERSTART (pymepix.SPIDR.spidrcmds.SpidrCmds attribute), 31 | at- | CMD_RESET_COUNTERS (pymepix.SPIDR.spidrcmds.SpidrCmds attribute), 31 |
| CMD_GET_SLVSCONFIG (pymepix.SPIDR.spidrcmds.SpidrCmds attribute), 31 | at- | CMD_RESET_DEVICES (pymepix.SPIDR.spidrcmds.SpidrCmds attribute), 31 |
| CMD_GET_SOFTWVERSION (pymepix.SPIDR.spidrcmds.SpidrCmds attribute), 31 | at- | CMD_RESET_MODULE (pymepix.SPIDR.spidrcmds.SpidrCmds attribute), 31 |
| CMD_GET_SPIDR_ADC (pymepix.SPIDR.spidrcmds.SpidrCmds attribute), 31 | at- | CMD_RESET_PIXELS (pymepix.SPIDR.spidrcmds.SpidrCmds attribute), 31 |
| CMD_GET_SPIDRREG (pymepix.SPIDR.spidrcmds.SpidrCmds attribute), 31 | at- | CMD_RESET_TIMER (pymepix.SPIDR.spidrcmds.SpidrCmds attribute), 31 |
| CMD_GET_STARTOPTS (pymepix.SPIDR.spidrcmds.SpidrCmds attribute), 31 | at- | CMD_RESTART_TIMERS (pymepix.SPIDR.spidrcmds.SpidrCmds attribute), 31 |
| CMD_GET_TIMER (pymepix.SPIDR.spidrcmds.SpidrCmds attribute), 31 | at- | CMD_SELECT_CHIPBOARD (pymepix.SPIDR.spidrcmds.SpidrCmds attribute), 31 |
| CMD_GET_TPNUMBER (pymepix.SPIDR.spidrcmds.SpidrCmds attribute), 31 | at- | CMD_SEQ_READOUT (pymepix.SPIDR.spidrcmds.SpidrCmds attribute), 31 |
| CMD_GET_TPPERIODPHASE (pymepix.SPIDR.spidrcmds.SpidrCmds attribute), 31 | at- | CMD_SET_BIAS_ADJUST (pymepix.SPIDR.spidrcmds.SpidrCmds attribute), 31 |
| CMD_GET_TRIGCONFIG (pymepix.SPIDR.spidrcmds.SpidrCmds attribute), 31 | at- | CMD_SET_BOARDDID (pymepix.SPIDR.spidrcmds.SpidrCmds attribute), 32 |
| | | CMD_SET_BUSY (pymepix.SPIDR.spidrcmds.SpidrCmds attribute), 32 |

| | | | |
|--|-----|---|-----|
| CMD_SET_CHIPBOARDID (<i>pymepix.SPIDR.spidrcmds.SpidrCmds</i> attribute), 32 | at- | CMD_SET_SERVERPORT (<i>pymepix.SPIDR.spidrcmds.SpidrCmds</i> attribute), 32 | at- |
| CMD_SET_CTPR (<i>pymepix.SPIDR.spidrcmds.SpidrCmds</i> attribute), 32 | | CMD_SET_SLVSCONFIG (<i>pymepix.SPIDR.spidrcmds.SpidrCmds</i> attribute), 32 | at- |
| CMD_SET_CTPR_LEON (<i>pymepix.SPIDR.spidrcmds.SpidrCmds</i> attribute), 32 | at- | CMD_SET_SPIDRREG (<i>pymepix.SPIDR.spidrcmds.SpidrCmds</i> attribute), 32 | at- |
| CMD_SET_DAC (<i>pymepix.SPIDR.spidrcmds.SpidrCmds</i> attribute), 32 | | CMD_SET_TIMEOFDAY (<i>pymepix.SPIDR.spidrcmds.SpidrCmds</i> attribute), 32 | at- |
| CMD_SET_DACS_DFLT (<i>pymepix.SPIDR.spidrcmds.SpidrCmds</i> attribute), 32 | at- | CMD_SET_TIMER (<i>pymepix.SPIDR.spidrcmds.SpidrCmds</i> attribute), 32 | at- |
| CMD_SET_EXTDAC (<i>pymepix.SPIDR.spidrcmds.SpidrCmds</i> attribute), 32 | | CMD_SET_TPNUMBER (<i>pymepix.SPIDR.spidrcmds.SpidrCmds</i> attribute), 32 | |
| CMD_SET_FANSPEED (<i>pymepix.SPIDR.spidrcmds.SpidrCmds</i> attribute), 32 | | CMD_SET_TPPERIODPHASE (<i>pymepix.SPIDR.spidrcmds.SpidrCmds</i> attribute), 32 | at- |
| CMD_SET_GENCONFIG (<i>pymepix.SPIDR.spidrcmds.SpidrCmds</i> attribute), 32 | at- | CMD_SET_TRIGCONFIG (<i>pymepix.SPIDR.spidrcmds.SpidrCmds</i> attribute), 32 | at- |
| CMD_SET_GPIO (<i>pymepix.SPIDR.spidrcmds.SpidrCmds</i> attribute), 32 | | CMD_STORE_ADDRPORTS (<i>pymepix.SPIDR.spidrcmds.SpidrCmds</i> attribute), 32 | |
| CMD_SET_GPIO_PIN (<i>pymepix.SPIDR.spidrcmds.SpidrCmds</i> attribute), 32 | | CMD_STORE_DACS (<i>pymepix.SPIDR.spidrcmds.SpidrCmds</i> attribute), 32 | at- |
| CMD_SET_HEADERFILTER (<i>pymepix.SPIDR.spidrcmds.SpidrCmds</i> attribute), 32 | at- | CMD_STORE_PIXCONF (<i>pymepix.SPIDR.spidrcmds.SpidrCmds</i> attribute), 32 | |
| CMD_SET_IPADDR_DEST (<i>pymepix.SPIDR.spidrcmds.SpidrCmds</i> attribute), 32 | at- | CMD_STORE_REGISTERS (<i>pymepix.SPIDR.spidrcmds.SpidrCmds</i> attribute), 32 | at- |
| CMD_SET_IPADDR_SRC (<i>pymepix.SPIDR.spidrcmds.SpidrCmds</i> attribute), 32 | at- | CMD_STORE_STARTOPTS (<i>pymepix.SPIDR.spidrcmds.SpidrCmds</i> attribute), 32 | |
| CMD_SET_LOGLEVEL (<i>pymepix.SPIDR.spidrcmds.SpidrCmds</i> attribute), 32 | | CMD_T0_SYNC (<i>pymepix.SPIDR.spidrcmds.SpidrCmds</i> attribute), 33 | at- |
| CMD_SET_OUTBLOCKCONFIG (<i>pymepix.SPIDR.spidrcmds.SpidrCmds</i> attribute), 32 | at- | CMD_TPX_POWER_ENA (<i>pymepix.SPIDR.spidrcmds.SpidrCmds</i> attribute), 33 | |
| CMD_SET_OUTPUTMASK (<i>pymepix.SPIDR.spidrcmds.SpidrCmds</i> attribute), 32 | at- | CMD_UPLOAD_PACKET (<i>pymepix.SPIDR.spidrcmds.SpidrCmds</i> attribute), 33 | at- |
| CMD_SET_PIXCONF (<i>pymepix.SPIDR.spidrcmds.SpidrCmds</i> attribute), 32 | | CMD_VALID_ADDRPORTS (<i>pymepix.SPIDR.spidrcmds.SpidrCmds</i> attribute), 33 | at- |
| CMD_SET_PLLCONFIG (<i>pymepix.SPIDR.spidrcmds.SpidrCmds</i> attribute), 32 | at- | CMD_VALID_DACS (<i>pymepix.SPIDR.spidrcmds.SpidrCmds</i> attribute), 33 | at- |
| CMD_SET_PWRPULSECONFIG (<i>pymepix.SPIDR.spidrcmds.SpidrCmds</i> attribute), 32 | at- | CMD_VALID_PIXCONF (<i>pymepix.SPIDR.spidrcmds.SpidrCmds</i> attribute), 33 | |
| CMD_SET_READOUTSPEED (<i>pymepix.SPIDR.spidrcmds.SpidrCmds</i> attribute), 32 | at- | CMD_VALID_REGISTERS (<i>pymepix.SPIDR.spidrcmds.SpidrCmds</i> attribute), 33 | at- |
| CMD_SET_SENSEDAC (<i>pymepix.SPIDR.spidrcmds.SpidrCmds</i> attribute), 32 | | | |

CMD_WRITE_FLASH (*pymepix.SPIDR.spidrcmds.SpidrCmd*
attribute), 33

columnTestPulseRegister
(*pymepix.SPIDR.spidrdevice.SpidrDevice*
attribute), 41

compute_timewalk () (in module *pymepix.util.timewalk*), 59

compute_timewalk_lookup () (in module *pymepix.util.timewalk*), 59

config (*pymepix.timepixdevice.TimepixDevice* attribute), 64

ConfigClassException, 63

configureStage () (*pymepix.processing.baseacquisition.AcquisitionStage*s (*pymepix.SPIDR.spidrcontroller.SPIDRController* method)), 51

connect_device () (*pymepix.processing.usbtrainid.USBTrainID*
method), 57

connect_timepix () (in module *pymepix.main*), 59

convertHtonl () (*pymepix.SPIDR.spidrcontroller.SPIDRController*
method), 35

convertNtohl () (*pymepix.SPIDR.spidrcontroller.SPIDRController*
method), 35

correct_global_time ()
(*pymepix.processing.logic.packet_processor.PacketProcessor*
method), 49

CpuToTpx (*pymepix.SPIDR.spidrcontroller.SPIDRController*
attribute), 33

create_socket_connection ()
(*pymepix.processing.udpsampler.UdpSampler*
method), 56

createConnection ()
(*pymepix.util.tcpsampler.TcpSampler* method), 58

cs_max_dist_tof (*pymepix.processing.logic.centroid_calculator.CentroidCalculator*
attribute), 46

cs_min_cluster_size
(*pymepix.processing.logic.centroid_calculator.CentroidCalculator*
attribute), 46

cs_sensor_size (*pymepix.processing.logic.centroid_calculator.CentroidCalculator*
attribute), 46

cs_tot_offset (*pymepix.processing.logic.centroid_calculator.CentroidCalculator*
attribute), 46

CycleOverflow (*pymepix.timepixdef.TimerOverflow*
attribute), 63

D

dacCodes () (*pymepix.config.defaultconfig.DefaultConfig*
method), 43

dacCodes () (*pymepix.config.sophyconfig.SophyConfig*
method), 44

dacCodes () (*pymepix.config.timepixconfig.TimepixConfig*
method), 45

DacRegisterCodes (class in *pymepix.timepixdef*), 61

data_thread () (*pymepix.pymepix_connection.PymepixConnection*
method), 61

dataCallback (*pymepix.pymepix_connection.PymepixConnection*
attribute), 60

datadrivenReadout ()
(*pymepix.SPIDR.spidrcontroller.SPIDRController*
method), 35

dbscan_clustering
(*pymepix.processing.logic.centroid_calculator.CentroidCalculator*
attribute), 47

Default (*pymepix.SPIDR.spidrdefs.SpidrReadoutSpeed*
attribute), 40

DefaultConfig (class in *pymepix.config.defaultconfig*), 43

DeviceCount (*pymepix.SPIDR.spidrcontroller.SPIDRController*
attribute), 36

deviceId (*pymepix.SPIDR.spidrdevice.SpidrDevice* attribute), 41

deviceIds (*pymepix.SPIDR.spidrcontroller.SPIDRController*
attribute), 36

deviceName (*pymepix.timepixdevice.TimepixDevice*
attribute), 65

DevicePort (*pymepix.SPIDR.spidrdevice.SpidrDevice*
attribute), 41

Disable (*pymepix.timepixdef.GrayCounter* attribute), 62

Disable (*pymepix.timepixdef.SuperPixel* attribute), 62

Disable (*pymepix.timepixdef.TestPulse* attribute), 63

disableExternalRefClock ()
(*pymepix.SPIDR.spidrcontroller.SPIDRController*
attribute), 56

disablePeriphClk80Mhz ()
(*pymepix.SPIDR.spidrcontroller.SPIDRController*
attribute), 36

DiscriminatorDigital

EmptyCalibpixdef.TestPulseDigAnalog
attribute), 63

enable (*pymepix.processing.basepipeline.BasePipelineObject*
attribute), 52

enable (pymepix.spidrcontroller.SPIDRController
attribute), 36

dvdd (*pymepix.SPIDR.spidrcontroller.SPIDRController*
attribute), 36

dvddNow (*pymepix.SPIDR.spidrcontroller.SPIDRController*
attribute), 36

E

enable (pymepix.spidrcontroller.SPIDRController
attribute), 36

Enable (*pymepix.timepixdef.GrayCounter* attribute), 62

Enable (*pymepix.timepixdef.SuperPixel* attribute), 62

Enable (*pymepix.timepixdef.TestPulse* attribute), 63

```

enable (pymepix.util.tcpsampler.TcpSampler attribute),           attribute), 41
      58                                         ExternalFallingTimer
enableDecoders () (pymepix.SPIDR.spidrcontroller.SPIDRController attribute), 41
      36                                         pymepix.SPIDR.spidrdefs.SpdrShutterMode attribute, 41
enableExternalRefClock ()                                         ExternalRisingFalling
      (pymepix.SPIDR.spidrcontroller.SPIDRController attribute), 41
enablePeriphClk80Mhz ()                                         ExternalRisingTimer
      (pymepix.SPIDR.spidrcontroller.SPIDRController attribute), 41
enablePolling () (pymepix.pymepix_connection.PymePixConnection attribute), 41
      61                                         pymepix.SPIDR.spidrcontroller.SPIDRController attribute
epsilon (pymepix.processing.logic.centroid_calculator.CentroidCalculator attribute), 37
      37                                         pymepix.SPIDR.spidrdefs.SpdrShutterCounter attribute
ERR_ADC_HARDW (pymepix.SPIDR.error.SPIDRErrorDefs attribute), 29
      29                                         filename (pymepix.config.sophyconfig.SophyConfig attribute), 44
ERR_DAC_HARDW (pymepix.SPIDR.error.SPIDRErrorDefs attribute), 29
      29                                         find_events_fast ()
ERR_FLASH_STORAGE                                         find_events_fast_post ()
      (pymepix.SPIDR.error.SPIDRErrorDefs attribute), 49
ERR_ILLEGAL_PAR (pymepix.SPIDR.error.SPIDRErrorDefs attribute), 29
      29                                         firmwareVersion (pymepix.SPIDR.spidrcontroller.SPIDRController attribute), 37
ERR_MON_HARDW (pymepix.SPIDR.error.SPIDRErrorDefs attribute), 29
      29                                         fpgaTemperature (pymepix.SPIDR.spidrcontroller.SPIDRController attribute), 37
ERR_MONITOR (pymepix.SPIDR.error.SPIDRErrorDefs attribute), 29
      29                                         frontEndAnalog (pymepix.timepixdef.TestPulseDigAnalog attribute), 63
ERR_MSG_LENGTH (pymepix.SPIDR.error.SPIDRErrorDefs attribute), 29
      29                                         getAdc () (pymepix.SPIDR.spidrcontroller.SPIDRController method), 58
ERR_NONE (pymepix.SPIDR.error.SPIDRErrorDefs attribute), 29
      29                                         getBuffers () (pymepix.processing.logic.packet_processor.PacketProcessor method), 49
ERR_NOT_IMPLEMENTED                                         getDac () (pymepix.SPIDR.spidrdevice.SpdrDevice method), 37
      (pymepix.SPIDR.error.SPIDRErrorDefs attribute), 29
ERR_SEQUENCE (pymepix.SPIDR.error.SPIDRErrorDefs attribute), 29
      29                                         getDacOut () (pymepix.SPIDR.spidrdevice.SpdrDevice method), 41
ERR_STR (pymepix.SPIDR.error.PymePixException attribute), 29
      29                                         getDevice () (pymepix.pymepix_connection.PymePixConnection method), 61
ERR_TPX3_HARDW (pymepix.SPIDR.error.SPIDRErrorDefs attribute), 29
      29                                         getLogger () (pymepix.core.log.Logger method), 45
ERR_UNKNOWN_CMD (pymepix.SPIDR.error.SPIDRErrorDefs attribute), 29
      29                                         getLogger () (pymepix.core.log.ProcessLogger method), 45
errorMessage () (pymepix.SPIDR.error.PymePixException method), 29
      29                                         getPixelConfig () (pymepix.SPIDR.spidrdevice.SpdrDevice method), 42
event_window (pymepix.processing.logic.packet_processor.PacketProcessor attribute), 49
      49                                         pymepix.util.tcpsampler.TcpSampler method), 58
EventData (pymepix.processing.datatypes.MessageType attribute), 53
      53                                         genConfig (pymepix.SPIDR.spidrdevice.SpdrDevice attribute), 41
EventiTot (pymepix.timepixdef.OperationMode attribute), 62
      62                                         get_useful_packets ()
External (pymepix.timepixdef.TestPulseGenerator attribute), 63
      63                                         (pymepix.processing.udpsampler.UdpSampler method), 56
ExternalFallingRising                                         get_useful_packets ()
      (pymepix.SPIDR.spidrdefs.SpdrShutterMode attribute)

```

getSpidrReg () (pymepix.SPIDR.spidrcontroller.SPIDRCController.krum (pymepix.timepixdevice.TimepixDevice method), 37
getStage () (pymepix.processing.baseacquisition.AcquisitionPipelinePixelDAC (pymepix.timepixdef.DacRegisterCodes attribute), 61
GrayCounter (class in pymepix.timepixdef), 62
grayCounter (pymepix.timepixdevice.TimepixDevice attribute), 65
gui_bar_fun (pymepix.post_processing.ProgressBar attribute), 59

H

handle () (pymepix.util.spidrDummyTCP.TPX3Handler method), 57
handle_events (pymepix.processing.logic.packet_processor attribute), 49
handle_lsb_time () (pymepix.processing.rawfilesampler.RawFileSampler method), 55
handle_msb_time () (pymepix.processing.rawfilesampler.RawFileSampler method), 55
handle_other () (pymepix.processing.rawfilesampler.RawFileSampler method), 55
hasOutput () (pymepix.processing.basepipeline.BasePipelineObject method), 64
headerFilter (pymepix.SPIDR.spidrdevice.SpidrDevice attribute), 42
HighSpeed (pymepix.SPIDR.spidrdefs.SpidrReadoutSpeed attribute), 40
humidity (pymepix.SPIDR.spidrcontroller.SPIDRController attribute), 37

I

Ibias_CP_PLL (pymepix.timepixdef.DacRegisterCodes attribute), 61
Ibias_DiscS1_OFF (pymepix.timepixdef.DacRegisterCodes attribute), 61
Ibias_DiscS1_OFF (pymepix.timepixdevice.TimepixDevice attribute), 64
Ibias_DiscS1_ON (pymepix.timepixdef.DacRegisterCodes attribute), 61
Ibias_DiscS1_ON (pymepix.timepixdevice.TimepixDevice attribute), 64
Ibias_DiscS2_OFF (pymepix.timepixdef.DacRegisterCodes attribute), 61
Ibias_DiscS2_OFF (pymepix.timepixdevice.TimepixDevice attribute), 64
Ibias_DiscS2_ON (pymepix.timepixdef.DacRegisterCodes attribute), 61
Ibias_DiscS2_ON (pymepix.timepixdevice.TimepixDevice attribute), 64
Ibias_Ikrum (pymepix.timepixdef.DacRegisterCodes attribute), 61

L

Ibias_PixelDAC (pymepix.timepixdevice.TimepixDevice attribute), 64
Ibias_Preamp_OFF (pymepix.timepixdef.DacRegisterCodes attribute), 61
Ibias_Preamp_OFF (pymepix.timepixdevice.TimepixDevice attribute), 64
Ibias_Preamp_ON (pymepix.timepixdef.DacRegisterCodes attribute), 61
Ibias_Preamp_ON (pymepix.timepixdevice.TimepixDevice attribute), 64
Ibias_TPbufferIn (pymepix.timepixdef.DacRegisterCodes attribute), 61
Ibias_TPbufferIn (pymepix.timepixdevice.TimepixDevice attribute), 64
Ibias_TPbufferOut (pymepix.timepixdef.DacRegisterCodes attribute), 61
Ibias_TPbufferOut (pymepix.timepixdevice.TimepixDevice attribute), 64
init_new_process ()
(pymepix.processing.pipeline_packet_processor.PipelinePacketProcessor method), 54
init_new_process ()
(pymepix.processing.rawfilesampler.RawFileSampler method), 55
init_new_process ()
(pymepix.processing.udpsampler.UdpSampler method), 56
init_new_process ()
(pymepix.util.tcp sampler.TcpSampler method), 58
Internal (pymepix.timepixdef.TestPulseGenerator attribute), 63
ipAddrDest (pymepix.SPIDR.spidrdevice.SpidrDevice attribute), 42
ipAddrSrc (pymepix.SPIDR.spidrdevice.SpidrDevice attribute), 42
isAcquiring (pymepix.pymepix_connection.PymepixConnection attribute), 61
isRunning (pymepix.processing.baseacquisition.AcquisitionPipeline attribute), 51
Left (pymepix.processing.logic.packet_processor.PixelOrientation attribute), 49
linkCounts (pymepix.SPIDR.spidrcontroller.SPIDRController attribute), 37
linkStatus (pymepix.SPIDR.spidrdevice.SpidrDevice attribute), 42

```

load_config()           (in      module
    pymepix.config.load_config), 43
loadConfig() (pymepix.timepixdevice.TimepixDevice
    method), 65
loadFile() (pymepix.config.sophyconfig.SophyConfig
    method), 44
localTemperature (pymepix.SPIDR.spidrcontroller.SPIDRController
    attribute), 37
Logger (class in pymepix.core.log), 45
LowSpeed (pymepix.SPIDR.spidrdefs.SpidrReadoutSpeed
    attribute), 40

```

M

```

main() (in module pymepix.config.sophyconfig), 44
main() (in module pymepix.main), 59
main() (in module pymepix.processing.baseacquisition),
    52
main() (in module pymepix.processing.basepipeline),
    53
main() (in module pymepix.processing.rawtodisk), 56
main() (in module pymepix.processing.udpsampler), 57
main() (in module pymepix.processing.usbtrainid), 57
main() (in module pymepix.SPIDR.spidrcontroller), 40
main() (in module pymepix.timepixdevice), 66
main() (in module pymepix.util.spidrDummyTCP), 57
main() (in module pymepix.util.spidrDummyUDP), 57
main() (in module pymepix.util.tcp sampler), 59
main_process() (in      module
    pymepix.processing.rawtodisk), 56

```

```

Mask (pymepix.timepixdef.GrayCounter attribute), 62
Mask (pymepix.timepixdef.OperationMode attribute), 62
Mask (pymepix.timepixdef.SuperPixel attribute), 62
Mask (pymepix.timepixdef.TestPulse attribute), 63
Mask (pymepix.timepixdef.TestPulseDigAnalog
    attribute), 63
Mask (pymepix.timepixdef.TestPulseGenerator
    attribute), 63
Mask (pymepix.timepixdef.TimeofArrivalClock
    attribute), 63
Mask (pymepix.timepixdef.TimerOverflow attribute), 63

```

```

maskPixels (pymepix.config.defaultconfig.DefaultConfig
    attribute), 43

```

```

maskPixels (pymepix.config.sophyconfig.SophyConfig
    attribute), 44

```

```

maskPixels (pymepix.config.timepixconfig.TimepixConfig
    attribute), 45

```

```

MessageType (class in pymepix.processing.datatypes),
    53

```

```

min_samples (pymepix.processing.logic.centroid_calculator.CentroidCalculator
    attribute), 47

```

```

MONITOR_ERR_STR (pymepix.SPIDR.error.PymePixException
    attribute), 29

```

N

```

Negative (pymepix.timepixdef.Polarity attribute), 62
numBlobProcesses (pymepix.processing.acquisition.CentroidPipeline
    attribute), 50
numDevices (pymepix.pymepix_connection.PymepixConnection
    attribute), 61
numToProcess (pymepix.processing.baseacquisition.AcquisitionStage
    attribute), 51

```

O

```

Open (pymepix.SPIDR.spidrdefs.SpidrShutterMode
    attribute), 41
open_file() (pymepix.processing.rawtodisk.Raw2Disk
    method), 55
open_output_file() (in      module
    pymepix.util.storage), 58
OpenFileCommand (pymepix.processing.datatypes.MessageType
    attribute), 53
openShutter () (pymepix.SPIDR.spidrcontroller.SPIDRController
    method), 37
OperationMode (class in pymepix.timepixdef), 62
operationMode (pymepix.timepixdevice.TimepixDevice
    attribute), 65
orientPixels () (pymepix.processing.logic.packet_processor.PacketPro
    method), 49
outBlockConfig (pymepix.SPIDR.spidrdevice.SpidrDevice
    attribute), 42
outfile_name (pymepix.processing.udpsampler.UdpSampler
    attribute), 56
outfile_name (pymepix.util.tcp sampler.TcpSampler
    attribute), 58
outputQueue (pymepix.processing.baseacquisition.AcquisitionStage
    attribute), 51
outputQueues (pymepix.processing.basepipeline.BasePipelineObject
    attribute), 53

```

P

```

PacketProcessor (class      in
    pymepix.processing.logic.packet_processor),
    48
PacketType (class in pymepix.timepixdef), 62
parseDAC () (pymepix.config.sophyconfig.SophyConfig
    method), 44
parsePixelConfig ()
    (pymepix.config.sophyconfig.SophyConfig
    method), 44
pauseHeartbeat () (pymepix.timepixdevice.TimepixDevice
    method), 65
pauseReadout (pymepix.SPIDR.spidrcontroller.SPIDRController
    method), 37
perform() (pymepix.clustering.cluster_stream.ClusterStream
    method), 43
perform_centroiding ()
    (pymepix.processing.logic.centroid_calculator.CentroidCalculator
    method), 43

```

method), 47
perform_centroiding_cluster_stream() *(pymepix.processing.logic.centroid_calculator.CentroidCalculator, 47)*
post_process () (pymepix.processing.logic.packet_processor.PacketProcessor, 49)
method), 49
method), 50
post_run () (pymepix.processing.basepipeline.BasePipelineObject, 53)
(pymepix.processing.logic.centroid_calculator.CentroidCalculator, 53)
method), 53
post_run () (pymepix.processing.pipeline_packet_processor.PipelinePacketProcessor, 54)
(pymepix.processing.logic.centroid_calculator.CentroidCalculator, 54)
method), 54
post_run () (pymepix.processing.rawfilesampler.RawFileSampler, 55)
(pymepix.processing.logic.centroid_calculator.CentroidCalculator, 55)
method), 55
post_run () (pymepix.processing.udpsampler.UdpSampler, 56)
(pymepix.processing.logic.centroid_calculator.CentroidCalculator, 56)
method), 56
post_run () (pymepix.util.tcp sampler.TcpSampler, 58)
PhaseShiftedGray (pymepix.timepixdef.TimeofArrivalClock, 58)
attribute), 63
powerPulseConfig (pymepix.SPIDR.spidrdevice.SpidrDevice, 58)
PipelineCentroidCalculator (class in pymepix.processing.pipeline_centroid_calculator, 42)
attribute), 42
pre_process () (pymepix.processing.logic.centroid_calculator.CentroidCalculator, 47)
54
PipelinePacketProcessor (class in pymepix.processing.pipeline_packet_processor, 49)
attribute), 49
pre_process () (pymepix.processing.logic.processing_step.ProcessingStep, 50)
method), 50
Pixel (pymepix.timepixdef.PacketType attribute), 62
PixelData (pymepix.processing.datatypes.MessageType, 53)
attribute), 53
pre_run () (pymepix.processing.basepipeline.BasePipelineObject, 53)
pixelMask (pymepix.timepixdevice.TimepixDevice, 65)
attribute), 65
pre_run () (pymepix.processing.pipeline_packet_processor.PipelinePacketProcessor, 54)
PixelOrientation (class in pymepix.processing.logic.packet_processor, 49)
attribute), 49
pre_run () (pymepix.processing.rawfilesampler.RawFileSampler, 55)
pixelPacketCounter (pymepix.SPIDR.spidrdevice.SpidrDevice, 56)
attribute), 42
pre_run () (pymepix.util.tcp sampler.TcpSampler, 58)
PixelPipeline (class in pymepix.processing.acquisition, 50)
attribute), 50
pre_pressure (pymepix.SPIDR.spidrcontroller.SPIDRController, 37)
pixelTest (pymepix.timepixdevice.TimepixDevice, 65)
attribute), 65
process () (pymepix.processing.basepipeline.BasePipelineObject, 53)
method), 53
pixelThreshold (pymepix.timepixdevice.TimepixDevice, 65)
attribute), 65
process () (pymepix.processing.logic.centroid_calculator.CentroidCalculator, 47)
PLL_Vcntrl (pymepix.timepixdef.DacRegisterCodes, 61)
attribute), 61
process () (pymepix.processing.logic.packet_processor.PacketProcessor, 49)
pllConfig (pymepix.SPIDR.spidrdevice.SpidrDevice, 42)
attribute), 42
process () (pymepix.processing.logic.processing_step.ProcessingStep, 50)
Polarity (class in pymepix.timepixdef, 62)
polarity (pymepix.timepixdevice.TimepixDevice, 62)
attribute), 65
process () (pymepix.processing.pipeline_centroid_calculator.PipelineCentroidCalculator, 54)
poll () (pymepix.pymepix_connection.PymepixConnection, 55)
method), 61
process () (pymepix.processing.pipeline_packet_processor.PipelinePacketProcessor, 55)
PollBufferEmpty, 60
process_pixels () (pymepix.processing.logic.packet_processor.PacketProcessor, 49)
pollBufferLength (pymepix.pymepix_connection.PymepixConnection, 61)
attribute), 61
trigger1 () (pymepix.processing.logic.packet_processor.PacketProcessor, 49)
Positive (pymepix.timepixdef.Polarity attribute), 62
post_process () (in module pymepix.main, 59)
process_trigger2 () (pymepix.processing.logic.packet_processor.PacketProcessor, 49)
post_process () (pymepix.processing.logic.centroid_calculator.CentroidCalculator, 49)
method), 49

processes (*pymepix.processing.baseacquisition.AcquisitionStage*.SPIDR.error (module), 29
attribute), 51
ProcessingStep (class in *pymepix.processing.logic.processing_step*), 49
ProcessLogger (class in *pymepix.core.log*), 45
ProgressBar (class in *pymepix.post_processing*), 59
PulseCounter (*pymepix.SPIDR.spidrdefs.SpidrShutterMode*.*pymepix.util* (module), 59
attribute), 41
push_data () (*pymepix.processing.rawfilesampler.RawFileSampler*.*util.spidrDummyTCP* (module), 57
method), 55
pushOutput () (*pymepix.processing.basepipeline.BasePipelineObject*.*util.tcp sampler* (module), 58
method), 53
pymepix (module), 66
pymepix.clustering (module), 43
pymepix.clustering.cluster_stream (module), 43
pymepix.config (module), 45
pymepix.config.defaultconfig (module), 43
pymepix.config.load_config (module), 43
pymepix.config.sophyconfig (module), 44
pymepix.config.timepixconfig (module), 44
pymepix.core (module), 45
pymepix.core.log (module), 45
pymepix.main (module), 59
pymepix.post_processing (module), 59
pymepix.processing (module), 57
pymepix.processing.acquisition (module), 50
pymepix.processing.baseacquisition (module), 51
pymepix.processing.basepipeline (module), 52
pymepix.processing.datatypes (module), 53
pymepix.processing.logic (module), 50
pymepix.processing.logic.centroid_calculator (module), 46
pymepix.processing.logic.packet_processor (module), 48
pymepix.processing.logic.processing_step (module), 49
pymepix.processing.logic.shared_processing_parameters (module), 50
pymepix.processing.pipeline_centroid_calculator (module), 54
pymepix.processing.pipeline_packet_processor (module), 54
pymepix.processing.rawfilesampler (module), 55
pymepix.processing.rawtodisk (module), 55
pymepix.processing.udpsampler (module), 56
pymepix.processing.usbtrainid (module), 57
pymepix.pymepix_connection (module), 60
pymepix.SPIDR (module), 43
pymepix.SPIDR.spidrcmds (module), 29
pymepix.SPIDR.spidrcontroller (module), 33
pymepix.SPIDR.spidrdefs (module), 40
pymepix.SPIDR.spidrdevice (module), 41
pymepix.timepixdef (module), 61
pymepix.timepixdevice (module), 63
pymepix.util.spidrDummyUDP (module), 57
pymepix.util.storage (module), 58
pymepix.util.tcp sampler (module), 58
pymepix.util.timewalk (module), 59
PymepixConnection (class in *pymepix.pymepix_connection*), 60
PymePixException, 29

R

Raw2Disk (class in *pymepix.processing.rawtodisk*), 55
RawData (*pymepix.processing.datatypes.MessageType* attribute), 53
RawFileSampler (class in *pymepix.processing.rawfilesampler*), 55
readoutSpeed (*pymepix.SPIDR.spidrdevice.SpidrDevice* attribute), 42
record (*pymepix.processing.udpsampler.UdpSampler* attribute), 56
record (*pymepix.util.tcp sampler.TcpSampler* attribute), 58
refreshPixels () (*pymepix.timepixdevice.TimepixDevice* method), 65
reinitDevice () (*pymepix.SPIDR.spidrdevice.SpidrDevice* method), 42
reinitDevices () (*pymepix.SPIDR.spidrcontroller.SPIDRController* method), 38
setTemperature (*pymepix.SPIDR.spidrcontroller.SPIDRController* attribute), 38
request () (*pymepix.SPIDR.spidrcontroller.SPIDRController* method), 38
requestGetBytes () (*pymepix.SPIDR.spidrcontroller.SPIDRController* method), 38
requestGetInt () (*pymepix.SPIDR.spidrcontroller.SPIDRController* method), 38
requestGetIntBytes () (*pymepix.SPIDR.spidrcontroller.SPIDRController* method), 38
requestGetInts () (*pymepix.SPIDR.spidrcontroller.SPIDRController* method), 38
requestSetInt () (*pymepix.SPIDR.spidrcontroller.SPIDRController* method), 38
requestSetIntBytes () (*pymepix.SPIDR.spidrcontroller.SPIDRController* method), 38

method), 38
`requestSetInts()` (*pymepix.SPIDR.spidrcontroller.SPIDRController*)
method), 38
`reset()` (*pymepix.SPIDR.spidrdevice.SpidrDevice*)
method), 42
`resetCounters()` (*pymepix.SPIDR.spidrcontroller.SPIDRController*)
method), 38
`resetDevices()` (*pymepix.SPIDR.spidrcontroller.SPIDRController*)
method), 38
`resetModule()` (*pymepix.SPIDR.spidrcontroller.SPIDRController*)
method), 38
`resetPacketCounters()` (*pymepix.SPIDR.spidrcontroller.SPIDRController*)
method), 39
`resetPixelConfig()` (*pymepix.SPIDR.spidrdevice.SpidrDevice*)
method), 42
`resetPixels()` (*pymepix.SPIDR.spidrdevice.SpidrDevice*)
method), 42
`resetPixels()` (*pymepix.timepixdevice.TimepixDevice*)
method), 65
`resetTimers()` (*pymepix.SPIDR.spidrcontroller.SPIDRController*)
pymepix.SPIDR.spidrdevice.SpidrDevice
method), 39
`restartTimers()` (*pymepix.SPIDR.spidrcontroller.SPIDRController*)
method), 39
`resumeHeartbeat()` (*pymepix.timepixdevice.TimepixDevice*)
method), 65
`Right` (*pymepix.processing.logic.packet_processor.PixelOrientation*)
attribute), 49
`run()` (*pymepix.processing.basepipeline.BasePipelineObject*)
method), 53
`run()` (*pymepix.processing.rawfilesampler.RawFileSampler*)
method), 55
`run()` (*pymepix.processing.udpsampler.UdpSampler*)
method), 57
`run()` (*pymepix.processing.usbtrainid.USBTrainID*)
method), 57
`run()` (*pymepix.util.tcpsampler.TcpSampler*)
method), 59
`run_post_processing()` (*in module*
pymepix.post_processing), 59
S
`saveMask()` (*pymepix.config.sophyconfig.SophyConfig*)
method), 44
`saveToHDF5()` (*pymepix.processing.rawfilesampler.RawFileSampler*)
method), 55
`sequentialReadout()` (*pymepix.SPIDR.spidrcontroller.SPIDRController*)
method), 39
`serverPort` (*pymepix.SPIDR.spidrdevice.SpidrDevice*)
attribute), 42
setArgs() (*pymepix.processing.baseacquisition.AcquisitionStage*)
setBiasSupplyEnable()
(pymepix.SPIDR.spidrcontroller.SPIDRController
method), 39
(pymepix.SPIDR.spidrcontroller.SPIDRController
method), 39
(pymepix.timepixdevice.TimepixDevice
method), 65
(pymepix.SPIDR.spidrdevice.SpidrDevice
method), 42
setDac() (*pymepix.timepixdevice.TimepixDevice*)
method), 65
setDacDefault() (*pymepix.SPIDR.spidrdevice.SpidrDevice*)
method), 42
setEthernetFilter()
(pymepix.timepixdevice.TimepixDevice
method), 65
setExternalDac() (*pymepix.SPIDR.spidrdevice.SpidrDevice*)
method), 42
setHeaderFilter()
setPixelMask() (*pymepix.SPIDR.spidrdevice.SpidrDevice*)
method), 42
setPixelTestBit()
(pymepix.SPIDR.spidrdevice.SpidrDevice
method), 42
setPowerPulseEnable()
(pymepix.SPIDR.spidrcontroller.SPIDRController)
method), 39
setSenseDac() (*pymepix.SPIDR.spidrdevice.SpidrDevice*)
method), 42
setShutterTriggerConfig()
(pymepix.SPIDR.spidrcontroller.SPIDRController)
method), 39
setSinglePixelMask()
(pymepix.SPIDR.spidrdevice.SpidrDevice)
method), 42
setSinglePixelTestBit()
(pymepix.SPIDR.spidrdevice.SpidrDevice)
method), 42
setSinglePixelThreshold()
(pymepix.SPIDR.spidrdevice.SpidrDevice)
method), 42
setSpidrReg() (*pymepix.SPIDR.spidrcontroller.SPIDRController*)
method), 39
setTpPeriodPhase()
(pymepix.SPIDR.spidrdevice.SpidrDevice)

| | | | |
|--|--|---|------------|
| <i>method)</i> , 42 | | | |
| setTpxPowerPulseEnable () | (<i>pymepix.SPIDR.spidrcontroller.SPIDRController</i> . <i>SPIDR_IPMUX_CONFIG_I</i>) | <i>(pymepix.SPIDR.spidrdefs.SpidrRegs</i> | <i>at-</i> |
| <i>method), 39</i> | | <i>tribute), 40</i> | |
| setupAcquisition () | (<i>pymepix.timepixdevice.TimepixDevice</i>) | (<i>pymepix.SPIDR.spidrdefs.SpidrRegs</i> | <i>at-</i> |
| <i>method), 66</i> | <i>SPIDR_PIXEL_FILTER_I</i> | <i>tribute), 40</i> | |
| setupDevice () (<i>pymepix.timepixdevice.TimepixDevice</i>) | | (<i>pymepix.SPIDR.spidrdefs.SpidrRegs</i> | <i>at-</i> |
| <i>method), 66</i> | <i>SPIDR_PIXEL_PKT COUNTER_I</i> | <i>tribute), 40</i> | |
| SharedProcessingParameter (class in <i>pymepix.SPIDR.spidrdefs.SpidrRegs</i>) | | | <i>at-</i> |
| <i>pymepix.processing.logic.shared_processing_parameter</i>), 40 | <i>SPIDR_PIXEL_PKT COUNTER_OLD_I</i> | | |
| 50 | | | |
| shutterCounter (<i>pymepix.SPIDR.spidrcontroller.SPIDRController</i>) | <i>pymepix.SPIDR.spidrdefs.SpidrRegs</i> | | <i>at-</i> |
| <i>attribute), 39</i> | <i>tribute), 40</i> | | |
| shutterEnd (<i>pymepix.SPIDR.spidrdevice.SpidrDevice</i>) | <i>SPIDR_SHUTTERTRIG_CNT_I</i> | | |
| <i>attribute), 42</i> | (<i>pymepix.SPIDR.spidrdefs.SpidrRegs</i> | | <i>at-</i> |
| shutterStart (<i>pymepix.SPIDR.spidrdevice.SpidrDevice</i>) | | <i>tribute), 40</i> | |
| <i>attribute), 42</i> | <i>SPIDR_SHUTTERTRIG_CTRL_I</i> | | |
| shutterTriggerConfig | | (<i>pymepix.SPIDR.spidrdefs.SpidrRegs</i> | <i>at-</i> |
| (<i>pymepix.SPIDR.spidrcontroller.SPIDRController</i>) | | <i>tribute), 40</i> | |
| <i>attribute), 39</i> | <i>SPIDR_SHUTTERTRIG_DELAY_I</i> | | |
| ShutterTriggerCount | | (<i>pymepix.SPIDR.spidrdefs.SpidrRegs</i> | <i>at-</i> |
| (<i>pymepix.SPIDR.spidrcontroller.SPIDRController</i>) | | <i>tribute), 40</i> | |
| <i>attribute), 34</i> | <i>SPIDR_SHUTTERTRIG_FREQ_I</i> | | |
| ShutterTriggerCtrl | | (<i>pymepix.SPIDR.spidrdefs.SpidrRegs</i> | <i>at-</i> |
| (<i>pymepix.SPIDR.spidrcontroller.SPIDRController</i>) | | <i>tribute), 40</i> | |
| <i>attribute), 34</i> | <i>SPIDR_SHUTTERTRIG_LENGTH_I</i> | | |
| ShutterTriggerDelay | | (<i>pymepix.SPIDR.spidrdefs.SpidrRegs</i> | <i>at-</i> |
| (<i>pymepix.SPIDR.spidrcontroller.SPIDRController</i>) | | <i>tribute), 41</i> | |
| <i>attribute), 34</i> | <i>SPIDR_TDC_TRIGGERCOUNTER_I</i> | | |
| ShutterTriggerFreq | | (<i>pymepix.SPIDR.spidrdefs.SpidrRegs</i> | <i>at-</i> |
| (<i>pymepix.SPIDR.spidrcontroller.SPIDRController</i>) | | <i>tribute), 41</i> | |
| <i>attribute), 34</i> | <i>SPIDR_UDP_PKT COUNTER_I</i> | | |
| ShutterTriggerLength | | (<i>pymepix.SPIDR.spidrdefs.SpidrRegs</i> | <i>at-</i> |
| (<i>pymepix.SPIDR.spidrcontroller.SPIDRController</i>) | | <i>tribute), 41</i> | |
| <i>attribute), 34</i> | <i>SPIDR_UDP MON _PKTCOUNTER_I</i> | | |
| ShutterTriggerMode | | (<i>pymepix.SPIDR.spidrdefs.SpidrRegs</i> | <i>at-</i> |
| (<i>pymepix.SPIDR.spidrcontroller.SPIDRController</i>) | | <i>tribute), 41</i> | |
| <i>attribute), 35</i> | <i>SPIDR_UDP PAUSE _PKTCOUNTER_I</i> | | |
| slaveConfig (<i>pymepix.SPIDR.spidrdevice.SpidrDevice</i>) | | (<i>pymepix.SPIDR.spidrdefs.SpidrRegs</i> | <i>at-</i> |
| <i>attribute), 42</i> | <i>tribute), 41</i> | | |
| softwareVersion (<i>pymepix.SPIDR.spidrcontroller.SPIDRController</i>) | (class in <i>pymepix.SPIDR.spidrcmds</i>), 29 | | |
| <i>attribute), 39</i> | <i>SPIDRController</i> (class in <i>pymepix.SPIDR.spidrcontroller</i>), 33 | | |
| SophyConfig (class in <i>pymepix.config.sophyconfig</i>), 44 | <i>SpidrDevice</i> (class in <i>pymepix.SPIDR.spidrdevice</i>), 41 | | |
| SPIDR_CPU2TPX_WR_I | | | |
| (<i>pymepix.SPIDR.spidrdefs.SpidrRegs</i>) | | <i>SPIDRErrorDfs</i> (class in <i>pymepix.SPIDR.error</i>), 29 | |
| <i>attribute), 40</i> | | <i>spidrFanSpeed</i> (<i>pymepix.SPIDR.spidrcontroller.SPIDRController</i>) | |
| SPIDR_DEVICES_AND_PORTS_I | | <i>attribute), 40</i> | |
| (<i>pymepix.SPIDR.spidrdefs.SpidrRegs</i>) | | <i>SpidrReadoutSpeed</i> (class in <i>pymepix.SPIDR.spidrdefs</i>), 40 | |
| <i>attribute), 40</i> | | <i>SpidrShutterMode</i> (class in <i>pymepix.SPIDR.spidrdefs</i>), 41 | |
| SPIDR_ERR_STR (<i>pymepix.SPIDR.error.PymePixException</i>) | | | |
| <i>attribute), 29</i> | | | |
| SPIDR_FE_GTX_CTRL_STAT_I | | | |

```

stage (pymepix.processing.baseacquisition.AcquisitionStage      method), 42
      attribute), 51
stages (pymepix.processing.baseacquisition.AcquisitionPipeline triggerCounter
      attribute), 51
      (pymepix.SPIDR.spidrcontroller.SPIDRController
start () (pymepix.processing.baseacquisition.AcquisitionPipeline attribute), 35
      method), 51
testPixels (pymepix.config.defaultconfig.DefaultConfig
start () (pymepix.processing.baseacquisition.AcquisitionStage   attribute), 43
      method), 52
testPixels (pymepix.config.sophyconfig.SophyConfig
start () (pymepix.pymepix_connection.PymepixConnection   attribute), 44
      method), 61
testPixels (pymepix.config.timepixconfig.TimepixConfig
start () (pymepix.timepixdevice.TimepixDevice    attribute), 45
      method), 66
testPulse (class in pymepix.timepixdef), 62
start_recording () testPulse (pymepix.timepixdevice.TimepixDevice attribute), 66
      (pymepix.timepixdevice.TimepixDevice
      method), 66
TestPulseDigAnalog (class in pymepix.timepixdef), 63
startAutoTrigger () testPulseDigitalAnalog
      (pymepix.SPIDR.spidrcontroller.SPIDRController
      method), 40
startTrainID () (pymepix.processing.baseacquisition.AcquisitionStage, 66
      method), 52
      TestPulseGenerator (class in pymepix.timepixdef),
stop () (pymepix.processing.baseacquisition.AcquisitionPipeline 63
      method), 51
      testPulseGeneratorSource
stop () (pymepix.processing.baseacquisition.AcquisitionStage  (pymepix.timepixdevice.TimepixDevice
      method), 52
      attribute), 66
stop () (pymepix.pymepix_connection.PymepixConnection thresholdPixels (pymepix.config.defaultconfig.DefaultConfig
      method), 61
      attribute), 43
stop () (pymepix.timepixdevice.TimepixDevice thresholdPixels (pymepix.config.sophyconfig.SophyConfig
      method), 66
      attribute), 44
stop_recording () (pymepix.timepixdevice.TimepixDevice thresholdPixels (pymepix.config.timepixconfig.TimepixConfig
      method), 66
      attribute), 45
stopAutoTrigger () TimeofArrivalClock (class in pymepix.timepixdef),
      (pymepix.SPIDR.spidrcontroller.SPIDRController
      method), 40
      timeOfArrivalClock
StopOverflow (pymepix.timepixdef.TimerOverflow      (pymepix.timepixdevice.TimepixDevice
      attribute), 63
      attribute), 66
stopRaw2Disk () (pymepix.util.tcp sampler.TcpSampler TimepixConfig (class
      method), 59
      pymepix.config.timepixconfig), 44
stopTrainID () (pymepix.processing.baseacquisition.AcquisitionStage device (class in pymepix.timepixdevice), 63
      method), 52
      timer (pymepix.SPIDR.spidrdevice.SpidrDevice
store_centroid () (in module pymepix.util.storage), 42
      attribute), 58
      TimerOverflow (class in pymepix.timepixdef), 63
STORE_ERR_STR (pymepix.SPIDR.error.PymePixException timerOverflowControl
      attribute), 29
      (pymepix.timepixdevice.TimepixDevice
store_raw () (in module pymepix.util.storage), 58
      attribute), 66
store_toa () (in module pymepix.util.storage), 58
store_tof () (in module pymepix.util.storage), 58
SuperPixel (class in pymepix.timepixdef), 62
superPixel (pymepix.timepixdevice.TimepixDevice
      attribute), 66
SystemClock (pymepix.timepixdef.TimeofArrivalClock tpNumber (pymepix.SPIDR.spidrdevice.SpidrDevice
      attribute), 63
      attribute), 42
TpPeriodPhase (pymepix.SPIDR.spidrdevice.SpidrDevice
      attribute), 41
t0Sync () (pymepix.SPIDR.spidrdevice.SpidrDevice
      attribute), 58
TPX3_ERR_STR (pymepix.SPIDR.error.PymePixException
      attribute), 58

```

T

t0Sync () (pymepix.SPIDR.spidrdevice.SpidrDevice

attribute), 29
TPX3Handler (class in *pymepix.util.spidrDummyTCP*), 57
Trigger (*pymepix.timepixdef.PacketType* attribute), 62
TriggerData (*pymepix.processing.datatypes.MessageType* attribute), 54
triggers_processed (*pymepix.processing.logic.centroid_calculator.CentroidCalculator* attribute), 47

U

UdpMonPacketCounter (*pymepix.SPIDR.spidrcontroller.SPIDRController* attribute), 35
UdpPacketCounter (*pymepix.SPIDR.spidrcontroller.SPIDRController* attribute), 35
UdpPausePacketCounter (*pymepix.SPIDR.spidrcontroller.SPIDRController* attribute), 35
UdpSampler (class in *pymepix.processing.udpsampler*), 56
UnknownParameterTypeException, 50

W

Up (*pymepix.processing.logic.packet_processor.PixelOrientation* attribute), 49
update_timer () (*pymepix.timepixdevice.TimepixDevice* method), 66
update_to () (*pymepix.post_processing.ProgressBar* method), 59
updateBuffers () (*pymepix.processing.logic.packet_processor.PacketProcessor* method), 49
updateProgressBar () (in module *pymepix.post_processing*), 60
uploadPacket () (*pymepix.SPIDR.spidrdevice.SpidrDevice* method), 42
uploadPixelConfig () (*pymepix.SPIDR.spidrdevice.SpidrDevice* method), 42
uploadPixels () (*pymepix.timepixdevice.TimepixDevice* method), 66
USBTrainID (class in *pymepix.processing.usbtrainid*), 57

V

value (*pymepix.processing.logic.shared_processing_parameter.SharedProcessingParameter* attribute), 50
vdd (*pymepix.SPIDR.spidrcontroller.SPIDRController* attribute), 40
vddNow (*pymepix.SPIDR.spidrcontroller.SPIDRController* attribute), 40
Vfbk (*pymepix.timepixdef.DacRegisterCodes* attribute), 62
Vfbk (*pymepix.timepixdevice.TimepixDevice* attribute), 64

VPreamp_NCAS (*pymepix.timepixdef.DacRegisterCodes* attribute), 62
VPreamp_NCAS (*pymepix.timepixdevice.TimepixDevice* attribute), 64
Vthreshold_coarse (*pymepix.timepixdef.DacRegisterCodes* attribute), 62
Vthreshold_coarse (*pymepix.timepixdevice.TimepixDevice* attribute), 64
Vthreshold_fine (*pymepix.timepixdef.DacRegisterCodes* attribute), 62
Vthreshold_fine (*pymepix.timepixdevice.TimepixDevice* attribute), 64
VTP_coarse (*pymepix.timepixdevice.TimepixDevice* attribute), 64
VTP_fine (*pymepix.timepixdef.DacRegisterCodes* attribute), 62
VTP_fine (*pymepix.timepixdevice.TimepixDevice* attribute), 64